

45x Trinocular LED Stereo Optical Microscope Wide Field A23.1003

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A23.1003
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month




Product Specification

- Applications: Medical
- Certification: CE, Rohs
- Drawtube: Binocular, Trinocular
- Illumination: LED
- Product Name: Stereo Zoom Microscope
- Theory: Stereo Microscope
- Highlight: **stereo zoom microscopes,
stereo dissecting microscope**



More Images



Product Description

45x Trinocular LED Stereo Optical Microscope Wide Field A23.1003

Description:

- High quality, high function and excellent cost performance.
- Best choose with ergonomic design, sharp image display and wide viewing field.
- Providing a wide viewing field of $\Phi 20\text{mm}$ with the zoom ratio of 1:4.5.
- Ideal instrument in academic, medical and industrial field.

Specifications:

A23.1003 Stereo Microscope		
Head	Binocular Viewing Head, Inclined at 45°	
Eyepiece	Extra wide field eyepiece EW10x/20	
Zoom Objective	1x 4.5x	
Zoom Ratio	1:4.5	
Working Distance	97mm	
Illumination	LED Reflection illumination	
	LED Transmission illumination	
Optional Accessories		Item No.
Head	Binocular head, inclined at 60°	A23.1003-BH
	Trinocular head, inclined at 45°	A23.1003-TH
Eyepiece	WF15x/16	A51.1012-15
	WF20x/12	A51.1012-20
	10x Eyepiece with cross hair	A51.1014
Zoom Objective	0.7x 3x	A52.1021-0730
	2x	A52.1006-2
Auxiliary objective	0.75x	A52.1006-75
	0.63x	A52.1006-63
	0.5x	A52.1006-5
Photo Attachment		A55.1001
Video Attachment	Cmount 0.5x	A55.1002-05
	Cmount 1x	A55.1002-10


About us :

OPTO-EDU have been focusing in microscope and educational instrument export market for many years. Our customers are from U.S.A, U.K., Russia, Canada, Germany, Denmark, Poland, Sweden, U.A.E. Qatar, Saudia Arabia, Egypt, Argentina, Korea, Thailand, Brazil, and more countries. Concentrating on Optical and Educational products, we are committed to be clients' excellent consultant and supplier for products and services in this field. We strive to be a unique, long-term and excellent supplier for all kinds of microscope from China. We sincerely welcome customers from all over the world contact us at any time!



Trinocular with diopter adjustable

Stand with LED light



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China