



Industry Zoom Optical Stereo Microscope Wide Field Microscopes A23.2601

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A23.2601
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Industry
- Certification: CE, Rohs
- Drawtube: Binocular, Trinocular
- Illumination: LED
- Product Name: Stereo Zoom Microscope
- Theory: Stereo Microscope
- Highlight: **stereo zoom microscopes,
stereo dissecting microscope**



Product Description

Industry Zoom Optical Stereo Microscope Wide Field Microscopes A23.2601

Specifications:

A23.2601 Series Zoom Stereo Microscope				BL-A2	B-P3	BL-P4	T-A1	TL-A2
Eyepiece	WF10x/22mm, high eyepoint (with eyecup)							
Viewing Head & Zoom Body	Binocular head, inclined 45°, rotatable 360°, Interpupillary Distance 54mm-76mm, Zoom Ratio 1:6.6	●	●	●				
	Trinocular head, inclined 45°, rotatable 360°, Interpupillary Distance 54mm-76mm, Zoom Ratio 1:6.6					●	●	
Zoom range	0.7x-4.6x							
Objective	Built-in 1x							
Magnification	7x-46x							
Working Distance	110mm							
Focusing	Focus unit							
Stand	Pillar stand		●	●				
	Arm stand	●				●	●	
Illumination	Reflecting lamp 3W LED (optional : 12V/15W halogen)	●		●				●
	Transmitting lamp 3W LED (optional : 12V/15W halogen)	●		●				●
Stage	Black & white plate							
	Glass plate	●		●				●
A23.2601 Series Zoom Stereo Microscope Optical Data								
Auxiliary Objectives		0.5x	0.75x	1x(built in)		2x		
Working distance (mm)		189mm	122mm	110mm		36mm		
Eyepiece	Magnification	3.5 23x	5.3 34.5x	7 46x		14 92x		
WF10x/22mm	Field dia. (mm)	62.86 9.57	41.90 6.38	31.43 4.78		15.71 2.39		
Eyepiece	Magnification	5.3 34.5x	7.9 51.8x	10.5 69x		21 138x		
WF15x/17mm	Field dia.(mm)	48.57 7.39	32.38 4.93	24.29 3.70		12.14 1.85		
Eyepiece	Magnification	7 46x	10.5 69x	14 92x		28 184x		
WF20x/13mm	Field dia.(mm)	37.14 5.65	24.76 3.77	18.57 2.83		9.29 1.41		
A23.2601 Series Zoom Stereo Microscope Combination List								
				Item No.				
Eyepiece	WF10x/22mm			A51.2621-10				
Viewing Head & Zoom Body	Binocular head & zoom body (1:6.6) without eyepiece			A53.2611-2				
	Trinocular head & zoom body (1:6.6) without eyepiece			A53.2612-2				
Focus Unit	Focus unit without illumination			A54.2608-1				
	Focus unit with 3W LED reflecting illumination			A54.2608-2				
Stand	Arm stand with focus unit, without Light			A54.2609-01				
	Arm stand with focus unit and 3W LED illumination			A54.2609-02				
	Pillar stand without illumination			A54.2609-03				
	Pillar stand with 3W LED transmitting illumination			A54.2609-04				
	Stand with mirror			A54.2609-05				
	Universal stand			A54.2609-06				
	Rectangular big base,3W-LED transmitting illumination, base size: 330×300mm			A54.2609-07				
	Rectangular big base,without illumination, base size:330×300mm			A54.2609-08				
Complete set	is composed of Eyepiece + Head & zoom body + Focus unit + Stand							
A23.2601-BL-A2	A51.2621-10 + A53.2611-2 + A54.2608-2 + A54.2609-02							
A23.2601-B-P3	A51.2621-10 + A53.2611-2 + A54.2608-1 + A54.2609-03							
A23.2601-BL-P4	A51.2621-10 + A53.2611-2 + A54.2608-2 + A54.2609-04							
A23.2601-T-A1	A51.2621-10 + A53.2612-2 + A54.2608-1 + A54.2609-01							
A23.2601-TL-A2	A51.2621-10 + A53.2612-2 + A54.2608-2 + A54.2609-02							

And More Combinations As Per Your Request		
Optional Accessories		Item No.
Eyepiece	WF15x/17mm	A51.2621-15
	WF20x/13mm	A51.2621-20
Auxiliary Objective	0.5x W.D.189mm	A52.2621-05
	0.75x W.D. 122mm	A52.2621-75
	2x W.D. 36mm	A52.2621-20
CCD Adapter	0.35x, focusing adjustable	A55.2601-SZ35
	0.5x	A55.2601-SZ05
	0.7x	A55.2601-SZ07
	For NIKON or SONY or CANON or OLYMPUS SLR digital camera	A55.2621





A23.2601-B

A23.2601-BL



A23.2601-T



A23.2601-TL



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China