



1000X Digital Optical Microscope Halogen Lamp For Laboratory A31.0903

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A31.0903
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Laboratory
- Certification: CE, Rohs
- Drawtube: Binocular, Monocular, Trinocular
- Illumination: Halogen Lamp
- Max Magnification: 1000X
- Theory: Digital Optical Microscope
- Highlight: portable digital microscope,
optical microscope with digital camera



Product Description

Description:

- 3.0M Built-in Digital Camera, USB 2.0,
- Infinity Optical System
- PL10x/18mm, Wide Field, Wide Angle Eyepiece,
- Monocular, Binocular, Trinocular Head Optional
- Critical And Kohler Illumination Optional
- Drawer Type Lamp Holder For Lamp Replace Easily
- Dark Field, Polarizing, Phase Contrast All Optional

Specifications:

A31.0903 Digital Laboratory Biological Microscope		
Digital Camera	3.0M Built-in Digital Camera, USB 2.0, Automatic Exposure, Automatic White Balance, Display 90% Image From Field of View. Line By Line Scan.	
Head	Compensation Free Binocular Head, 30° Inclined,360° Rotatable,+/-5° Diopter Adjustable	
Eyepiece	Wide Field WF10x/18mm, High Eyepoint	
Nosepiece	Quadruple Nosepiece	
Objective	Infinity Plan Achromatic 4x/0.10, 10x/0.25, 40xs/0.65, 100xs/1.25(Oil)	
Stage	Double Layer Mechanical Stage, 132*140mm, Mvoing Range 50*76mm, Precision 0.1mm	
Focusing	Low Position Coaxial Coarse & Fine Focusing. Coarse Focusing Range 22mm, Fine Focusing Scale 0.002mm, With Tension Adjustment, With Up Limit	
Condenser	N.A.1.2/0.24 Swing Out Abbe Condenser, Iris Diaphargm, With Filter Holder	
Filter	Blue, Dia.32mm	
Light Source	Khohler Illumination, 6V/20W Halogen Lamp	
Power	90~240V Wide Range Voltage	
A31.0903 Digital Laboratory Biological Microscope Optional Accessories		
Eyepiece	WF10x/18mm with Reticle	A51.0905-1018R
	WF10x/18mm with Pointer	A51.0906-1018
	WF15x/13mm, High Eyepoint	A51.0902-1513
	WF16x/14mm, High Eyepoint	A51.0902-1614
	WF20x/10mm High Eyepoint	A51.0902-2010
Objective	Plan Achromatic 25xs/0.40	A52.0903-25
	Infinity Plan Achromatic 20xs/0.50	A52.0906-20
Adapter	0.5x C-Mount, For 1/2" CCD, Focus Adjustable	A55.0923-05
	0.35x C-Mount, For 1/3" CCD, Focus Adjustable	A55.0923-35
Polarizing	Polarizer Dia..32mm, Put on Collector	A5P.0916-P
	Analyzer Slide, Insert To Hole On Nosepiece	A5P.0916-A
	Rotatable Working Stage	A5P.0916-RW
Dark Field	Dark Field Condenser N.A.0.9, Dry	A5D.0911-D
	Dark Field Condenser N.A.1.25, Immersion	A5D.0911-W
Phase Contrast	Centering Telescope, Dia.23.2mm	A5C.0910-232
	5 Holes Turret Phase Contrast Condenser, N.A. 1.25	A5C.0912-PCC
	Phase Contrast Condenser for 10x	A5C.0911-10
	Phase Contrast Condenser for 20x	A5C.0911-20
	Phase Contrast Condenser for 40x	A5C.0911-40
	Phase Contrast Condenser for 100x	A5C.0911-100
	Phase Contrast Plan Achromatic Objective 10x	A5C.0931-10
	Phase Contrast Plan Achromatic Objective 20x	A5C.0931-20
	Phase Contrast Plan Achromatic Objective 40x	A5C.0931-40
	Phase Contrast Plan Achromatic Objective 100x	A5C.0931-100
Infinity Phase Contrast	Centering Telescope, Dia.23.2mm	A5C.0910-232
	5 Holes Turret Phase Contrast Condenser, N.A. 1.25	A5C.0912-PCCI
	Phase Contrast Condenser for 10x	A5C.0911-10
	Phase Contrast Condenser for 20x	A5C.0911-20
	Phase Contrast Condenser for 40x	A5C.0911-40
	Phase Contrast Condenser for 100x	A5C.0911-100
	Phase Contrast Infinity Plan Achromatic Objective 10x	A5C.0935-10
	Phase Contrast Infinity Plan Achromatic Objective 20x	A5C.0935-20
	Phase Contrast Infinity Plan Achromatic Objective 40x	A5C.0935-40
	Phase Contrast Infinity Plan Achromatic Objective 100x	A5C.0935-100



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China