



180X 540X USB Handheld Digital Microscope For Research A34.0601

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A34.0601
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Model Number: A34.4185
- Application: Research, Learning
- Certification: CE, Rohs
- Pixel: 1.3 Mega
- Theory: Handheld Digital Microscope
- Data Transfer: USB 2.0
- Highlight: **Portable Digital Microscope, hand held microscope**



More Images



Product Description

180X 540X USB Handheld Digital Microscope For Research A34.0601


Description:

- New Design in 2008!
- Magnification 30x - 540x.
- USB 2.0, including software
- Professional 1/2" industry grade 1.30M pixel CMOS image sensor

Specifications:

Specification	
Magnifications	30X,60X,180X,540X
Image Sensor	Professional 1/2" industry grade 1.3 mega pixels 1/2" Colorful CMOS
Pixel Size	5.2μm×5.2μm
Spectrum Response	400nm 1000nm
Sensitivity	1 Lux
Scanning Type	Progressive Scanning
Exposal Type	Electronic Rolling Shutter (ERS)
Resolution	1280×1024 pixels
Transmission Rate	16f/s@1280×1024 30f/s@640×480 30f/s@320×240
White Balance	Automatic/Manual Operations Supported
Exposal Control	Automatic/Manual Operations Supported
Image Output	USB2.0,480Mb/s
Power supply	USB 2.0 Power Supply, 800 mA
Dynamic Range	62dB
Programmable Control	Sizes, quality, luminance, gain, exposal time, contrast,saturation, resolution and Gamma of image
Illumination Control	Precise Laser locating, and top and bottom LED illumination grade control
Working Temperature	0°C~50°C
Software	Compatible with the operation systems of windows 2000, XP, VISTA and etc
Catalogue No.	Optional Accessories
A52.06SO-90	9x Auxillary Objective, magnification upto 540x, with working stage
A56.06LG-01	USB lateral LED illuminator







www.cnoec.com



www.cnoec.com



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China