



## Blue Portable 200x Hand Held Digital Microscope With White Light LED A34.4105

### Our Product Introduction

for more products please visit us on [cnoec.com](http://cnoec.com)

#### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A34.4105
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



#### Product Specification

- Model Number: A34.4105
- Application: Learning|Home
- Certification: CE|Rohs
- Magnification: 200x
- Theory: Handheld Digital Microscope
- Sensor: CMOS 0.3M
- Highlight: Portable Digital Microscope,  
Hand held Digital Microscope



#### More Images



## Product Description

Blue Portable 200x Hand Held Digital Microscope With White Light LED A34.4105

### Description:

**A34.4105 500x Wireless WIFI Digital Microscope** makes you have no limitation of distance; It can be used in a range within 30 to 80 meters. It is a perfect match with IOS mobile device (iPad, iPhone) and android mobile device (android tablet PC, android mobile phone). Please read below manual to get more detail about the microscope and follow the operation instructions.

- Textile industrial for textile inspection
- Printing inspection
- Industrial inspection: PCB, Precision machinery...
- Educational purpose
- Hair examination
- Skin examination
- Microbiological observation
- Jewelry & coin (Collections) inspection
- Visual Assistance

### Specifications:

A34.4105 500x Wireless WIFI Digital Microscope	
Magnification	Range 5X ~ 200X (Manually)
Light Source	Built-in 4~ 8 White-light LED and adjustable illumination ensure the magnified images are clear and bright
Focus Range	0mm ~ 40mm, other range can be customized
WIFI Standard	IEEE 802.11 b/g/n standards compliant
Resolution	Still Image Capture : 640*480, can be magnified and minified
Sensor	CMOS 0.3M
Controller	High Speed DSP (M-jpeg)
Take Images	Snap shot in build-in application program
Brightness	Manual adjustment
Power Supply	Rechargeable 1000mA self-protection lithium battery
Support O/S	WIFI capable, such as IOS and Android
Language	English (other language can be customized)
Application Software	CoolingTech
Illumination Scope	0 ~ 30000LUX
Product Dimension	142 mm ( L ) X 37 mm ( R )
Product Net Weight	100g





A34.4105 500x Wireless WIFI Digital Microscope



Power & WIFI Indicator, Micro USB Port



8 LEDs Light Source










A34.4105 500x Wireless WIFI Digital Microscope



**Opto-Edu (Beijing) Co., Ltd.**

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China