



## 1000x Multi Viewing Microscope 2 Position With Halogen Lamp A17.1101

Our Product Introduction

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE,
- Model Number: A17.1101
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Illumination: Halogen Lamp
- Product Name: Microscope
- Head: Binocular Head
- Eyepiece: WF10X/22mm
- Focusing: Coaxial Coarse & Fine Focusing
- Condenser: Abbe N.A. 1.25, With Iris Diaphragm And Filter
- Highlight: **dual head microscope, Multi View Microscope**



for more products please visit us on [cnoec.com](http://cnoec.com)

## Product Description

### 1000x Multi Viewing Microscope 2 Position With Halogen Lamp A17.1101

A17.1101 Multi-Viewing Microscope, 2 Position		A	B
Head	Dual Viewing Compensation Free Binocular Head, Interpupillary Distance 55-75mm		
Indicator	LED Cursor Pointer	-	o
Eyepiece	WF10X/22mm, 4 pcs		
Objective	Infinity Achromatic 4x, 10x, 40x(s), 100x(s, oil)		
Nosepiece	Quadruple Revolving Nosepiece		
Stage	Double Layer Mechanical Stage, Size 180x150mm, Moving Range 70x50mm		
Condenser	Abbe Condenser N.A.=1.25, with Iris Diaphragm and Filter		
Focusing	Coaxial Coarse & Fine Focusing Knobs, Scale:0.002mm		
Light Source	Halogen Lamp 6V 20W, Continue Intensity Control, Kohler Illumination		



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ [sale@optoedu.com](mailto:sale@optoedu.com)

🌐 [cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China