



WIFI 10X - 200X Handheld digital microscope For iPhone / iPad / PC / Android

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A34.5002
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Highlight: Hand held Digital Microscope, hand held microscope



More Images



Product Description

WIFI 10X - 200X Handheld digital microscope For iPhone / iPad / PC / Android

A34.5002 WIFI Miroscope, iPhone/iPad/Pc/Android, 10X-200X

- Magnification 10x~200x
- WIFI Transmission Rate 150Mbps
- Light Source 6 LEDs, Brightness Adjustable
- Support iPhone/iPad/Pc/Android/Mac Book
- It can work on Windows/Mac/iOS (iPhone, iPad)/Android (smart phone, Tablets)
- It has measurement function (by software) when working on Windows/Mac
- It comes with polarizer (which can reduce reflection when used to observe shiny objects like jewelry, metal etc)
- Professional Stand, very suitable for observation

A34.5002 WIFI Miroscope, iPhone/iPad/Pc/Android, 10X-200X	
Resolution	1.0M
Senser	1/4" CMOS, Color
Magnification	10x~200x
Video Resolution	720P (1280*720)
Video Format	H.264
Frame Rate	Wireless: 720P: 15fps, VGA: 25fps
Photo Resolution	1280*720
Snapshot	By Software
Picture Format	JPG
Wireless	WIFI Direct
Transmission Range	20M, Open Sight
Transmission Rate	150Mbps
Power Supply	DC5V, 1A, by Mini USB
Battery	Rechargeable Battery, Lithium 1300mAh, Can Be Charged By Adaptor/USB
Battery Working Time	Around 75 Minutes
Operation Temperature	0~65°
Light Source	6 LEDs, Brightness Adjustable
Polarizer	Yes
Operation System	iOS (5.1 or Later) Android Windows XP/Win7/Win8 Mac Book (10.7 or Later)
Measurement	Yes, For Windows and Mac Book
Software	MustCapture
Accessories	Professional Stand, Adapter, USB Cable, Softare CD, Calibration Charge
Dimension	145(L) x 36(D)mm
Weight	130g



A34.5002 WIFI Microscope, iPhone/iPad/Pc/Android, 10X-200X



www.cnoec.com





A34.5002 WIFI Microscope, With Stand Height Adjustable





A34.5002 WIFI Microscope, With Stand Height Adjustable



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China