



## 200x Wide Field Research Forensic Comparison Microscope A18.1801

### Our Product Introduction

#### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A18.1801
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



#### Product Specification

- Applications: Student, School, Research
- Certification: CE
- Theory: Comparison Microscope
- Highlight: **transmitted light microscope, binocular optical microscope**

## Product Description

### 200x Wide Field Research Forensic Comparison Microscope A18.1801

A18.1801 comparison microscope, is suitable to be used for forensic science education & school, to give an idea and training to the students what is a comparison microscope and learn how to use it. It is simple but including full comparison system function, and with steady quality at very competitive price!

- Total magnification 2.4x-200x
- 3 level change over objective 0.8x, 2x, 5x
- Air blast cooling light source provide high brightness
- Full left/right, split-image or overlapping image comparison
- Camera and video adapter included
- Magnification correcting system built-in

A18.1801 Comparison Microscope Specification	
Magnification	2.4x~200x
Eyepiece	WF10X/18mm
	WF20X/14mm
Objective	0.8x, 2x, 5x objective, with magnification correcting system
	Auxiliary lens 0.3x
	Auxiliary lens 2x
Bullet Holder	Small bullet holder set
Working Stage	Small working stage
Light Source	50W Air blast cooling light source provide high brightness
	Red/Green Filter
Adapter	SLR Camera Adapter
	CCD Adapter
Instrument Size	29mm*23mm*58mm

Optional Accessories		Item No.
Bullet Holder	Motorized bullet holder set	A54.1803
Light Source	150W cold light source, soft fiber Dia.12mm with positioning ruler	A56.1804
	Hard Optical Fiber, pair, for quick adjustment of light spot	A56.1805
	Transmit light lens	A56.1806
	Coaxial light lens	A56.1807
	Polarizing lens	A56.1808
	Ultraviolet light lens	A56.1810
Adapter	Digital camera adapter for Canon A650	A55.1803
Software	Professional forensic image analysis software <b>Crime Image</b>	A30.1801





0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China