



Zoom Stereo Optical Microscope With Tilting Head , 0.8 - 8x , A23.1005

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE,
- Model Number: A23.1005
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Industry, Medical
- Certification: CE
- Illumination: LED
- Theory: Stereo Microscope
- Highlight: **stereo dissecting microscope, stereo binocular microscopes**



More Images



Product Description

Zoom Stereo Optical Microscope With Tilting Head , 0.8 - 8x , A23.1005

A23.1005 Zoom Stereo Microscope With Tilting Head , the newest advanced zoom stereo microscope, can be used extensively for scientific research, industry, biological and medical field.

- Parallel Optical System, Module Design to Add Function Easily
- Tilting Head 0°-35° Viewing Angle Adjustable for User
- Large Zoom Ratio Up to 1:10
- Transmit & Reflect LED Illumination, Brightness Adjustable, Working Life 6000 Hours
- Complete Accessories Include Polarizing, Dark Field, Fluorescent Attachment
- Plan Apochromatic Objective 0.5x Optional

A23.1005 Zoom Stereo Microscope With Tilting Head , the newest advanced zoom stereo microscope, can be used extensively for scientific research, industry, biological and medical field.

- Parallel Optical System, Module Design to Add Function Easily
- Tilting Head 0°-35° Viewing Angle Adjustable for User
- Large Zoom Ratio Up to 1:10
- Transmit & Reflect LED Illumination, Brightness Adjustable, Working Life 6000 Hours
- Complete Accessories Include Polarizing, Dark Field, Fluorescent Attachment
- Plan Apochromatic Objective 0.5x Optional

A23.1005 Zoom Stereo Microscope, With Tilting head		A	B	C
Optical System	Parallel Optical System, Module Design to Add Function Easily			
Head	Tilting Head 0°-35° Viewing Angle Adjustable for User			
Eyepiece	Extra Wide Field Eyepiece EWF10x/22mm			
Zoom Body	0.8~5.0x, Zoom Ratio 1:6	o		
	0.8~6.4x, Zoom Ratio 1:8		o	
	0.8~8.0x, Zoom Ratio 1:10			o
Auxiliary Lens	Plan Achromatic Lens 1.0x			
Working Distance	78mm			
Focusing	Coarse Focusing Konb, Range 105mm	o	o	
	Coaxial Coarse & Fine Focusing Knob, Fine Division 0.015mm, Coarse Stroke 38.4mm/Rotation, Fine Stroke 1.5mm/Rotation, Moving Range 105mm			o
Illumination	Transmit & Reflect LED Illumination, Brightness Adjustable, Working Life 6000 Hours			
A23.1005 Zoom Stereo Microscope, With Tilting head		Cata. No.		
Eyepiece	WF15x/16mm	A51.1027-1516		
	WF20x/12mm	A51.1027-2012		
	WF30x	A51.1027-30		
Auxiliary Lens	Plan Apochromatic Lens 0.5x	A52.1023-05A		
	Plan Achromatic Lens 1.0x	A52.1023-10		
	Plan Achromatic Lens 2.0x	A52.1023-20		
	Achromatic Lens 0.3x	A52.1023-03		
	Achromatic Lens 0.5x	A52.1023-05		
Illumination	Coaxial Illumination Set: --Cold Light Power Supplier, --Dual Light Guide, --Coaxial Episcopic Illuminator, --Iris Diaphragm	A56.1021		
	Iris Diaphragm	A56.1021-2		
Working Stage	Temperature Control Stage	A54.1004		
	Movable Working Stage	A54.1007		
Stand Click to See	Universal Table Stand	A54.1011		
		A54.1012		
		A54.1013		
		A54.1014		
	Track Stand, Transmit & Reflect Light Source With Coarse & Fine Focusing Optional, Model -F, Like A54.1019-F	A54.1019		
		A54.1020		
		A54.1021		
Fluorescent	EPI-Fluorescent Attachment For Stereo Microscope	A5F.1020		
Dark Field	Dark Field Attachment For Stereo Microscope	A5D.1001		
Polarizing	Simple Polarizing Attachment For Stereo Microscope	A5P.1003		
Adapter	Photo Tube, Beamsplitter 1 Port	A55.1003-S1		
	Photo Tube, Beamsplitter 2 Port	A55.1003-S2		
	Photo Adapter, 2.0x, For Nikon or Canon Digital SLR Camera	A55.1020-2X		
	C-mount 0.5x	A55.1021-05		
	C-mount 1.0x	A55.1021-10		



A23.1005 Titling Head Zoom Stereo Microscope, 0.8x-5x, 6.4x, 8x



Titling Head 0°-35° Viewing Angle Adjustable



**Plan 1x Objective Lens
Up & Bottom Light Source, Brightness Adjustable**

**Parallel Optical System,
Module Design to Add Function Easily**





A23.1005 Standard Set

A23.1005 Standard Set With
A55.1003-S1 Photo Tube, Beamsplitter 1 Port





A23.1005 Zoom Stereo Microscope
Tilting Head 0°-35° Viewing Angle Adjustable



A23.1005 Zoom Stereo Microscope
Plan 1x Objective Lens
Up & Bottom Light Source, Brightness Adjustable

Parallel Optical System,
Module Design to Add Function Easily



A23.1005 Standard Set

**A23.1005 Standard Set With
A55.1003-S1 Photo Tube, Beamsplitter 1 Port**



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China