



20x - 40x Digital Stereo Microscope A32.1201 With Inclined Binocular Head

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A32.1201
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Certification: CE
- Drawtube: Binocular
- Illumination: Halogen Lamp
- Highlight: **portable digital microscope,**
optical microscope with digital camera

Product Description

20x - 40x Digital Stereo Microscope A32.1201 With Inclined Binocular Head


A32.1201 Digital Stereo Microscope		
Head	Inclined Binocular Head,360° Rotatable,Built-in CMOS(1.3M Pixel)	
Eyepiece	WF10X / Φ 20	
Objective	2x/4x	
Focusing	Coarse Focusing Adjustment	
Incident Light	Halogen Lamp 12V/10W,Brightness Adjustable	
Transmitted Light	Halogen Lamp 12V/10W,Brightness Adjustable	
Optional Accessories		Item No.
Eyepiece	WF5X/ Φ 16	A51.1201-05
	WF15X/ Φ 15	A51.1201-15
	WF20X/ Φ 10	A51.1201-10





 0086 13911110627

 sale@optoedu.com

 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China