



10x Stereo Optical Microscope Binocular WF10x A22.1201 With CE

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A22.1201
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Certification: CE
- Drawtube: Binocular
- Theory: Stereo Microscope
- Highlight: **stereo zoom microscopes,**
stereo binocular microscopes

for more products please visit us on cnoec.com

Product Description

10x Stereo Optical Microscope Binocular WF10x A22.1201 With CE

Standard Outfit:

1. Sharp stereo image over a wide field view
2. The left ocular-tube with diopter adjustment 5dp. The interpupillary distance is between 55-75mm
3. Dovetail inserted style objective with locked-in device ensures easy replacement
4. Power supply 110V/60Hz or 220V/50Hz with attached illumination AC 12V/10W or DC 6V/6W

Optional accessories:

1. Eyepiece WF5x(220mm), WF10x(20mm) WF15x(15mm), WF20x(10mm)
2. Objective 0.5x(250mm) 1x(230mm)
3. Incandescent bulb 12V/10W

A22.1201 Stereo Microscope									
Head	The Left Ocular-tube with diopter adjustment 5dp, the Interpupillary Distance:55-75mm								
Eyepiece	WF10X								
Objective	1X								
Working Distance	230mm								
Light Source	Incandescent Lamp 12V/10W 6V/6W								
Optional Accessories									Item No.
Eyepiece	WF5X								A51.1201-05
	WF15X								A51.1201-15
	WF20X								A51.1201-20
Objective	0.5X								A52.1281-05
Technical Data									
Objective	EYEPIECES								Working Distance
	5X		10X		15X		20X		
	Mag.	Field of View	Mag.	Field of View	Mag.	Field of View	Mag.	Field of View	
0.5X	2.5X	40.0mm	5X	40.0mm	7.5X	30.0mm	10X	20.0mm	250mm
1X	5X	20.0mm	10X	20.0mm	15X	15mm	20X	10.0mm	230mm



A22.1201-2D

A22.1201-6D



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China