



Binocular Head Plan Achromatic Objective Biological Compound Microscope 1000X A12.0201

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A12.0201
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Certification: CE
- Drawtube: Binocular, Trinocular
- Illumination: Halogen Lamp
- Max Magnification: 1000X
- Highlight: **optical compound microscope, wide field microscope**



Product Description

Binocular Head ,Plan Achromatic Objective Biological Compound Microscope 1000X A12.0201

This microscope equipped with plan achromatic objective and wide field eyepieces, can provide clear image by easy operation.

- With plan objectives
- Binocular or trinocular head inclined 30° available
- Coaxial coarse/fine focusing system with tension adjustable and up stop
- Kohler illumination with 6V 20W halogen lamp


Specification	A12.0201-A1	A12.0201-A2	A12.0201-B1	A12.0201-B2
Eyepiece	WF10X(Φ20mm)			
Objective	Plan 4X/0.10		Plan 4X/0.10	
	Plan 10X/0.25		Plan 10X/0.25	
			Plan 20X/0.40	
	Plan 40X/0.65(Spring)		Plan 40X/0.65(Spring)	
	Plan 100X/1.25(Spring, oil)		Plan 100X/1.25(Spring, oil)	
Head	Binocular	Trinocular	Binocular	Trinocular
Focus System	Coaxial coarse/fine focus system, Minimum division of fine focusing: 2um With tensional adjustable and up stop,			
Nosepiece	Quadruple backward ball bearing inner locating		Quintuple backward ball bearing inner locating	
Stage	Double layer mechanical, size 160x140mm, moving range 75x50mm			
Condenser	Abbe N.A.1.25 Rack & pinion adjustment			
Filter	Blue filter			
	Frosted filter			
Collector	For halogen lamp			
Light source	6V 20W halogen lamp with brightness control			

Optional Accessories		Item No.
Eyepiece	Wide field WF16x/11mm	A51.0203-16A
	Dividing 10x/20mm, 0.1mm/Div	A51.0205-10
Filter	Green filter	A56.0209-02
	Yellow filter	A56.0209-03
Condenser	Dark field, dry	A5D.0210-D
	Dark field, immersion	A5D.0210-I
Phase Contrast Attachment	Centering telescope	A5C.0212-1
	Objective, Phase contrast plan 10x/0.25 PHP	A5C.0230-P10
	Objective, Phase contrast plan 20x/0.40 PHP(Spring)	A5C.0230-P20
	Objective, Phase contrast plan 40x/0.65 PHP(Spring)	A5C.0230-P40
	Objective, Phase contrast plan 100x/1.25 PHP (Spring,oil)	A5C.0230-P100
	PH-I Turnplate phase contrast condenser	A5C.0210-1
	PH-II Turnplate phase contrast condenser	A5C.0210-4
	Flashboard phase contrast set	A5C.0210-2
	Pull board phase contrast set	A5C.0210-3
	Binocular head inclination of 30°	A53.0203-111
Infinity Optical System	Trinocular head inclination of 30°	A53.0208-1
	Centering telescope	A5C.0212-2
	Objective, infinity plan achromatic 4x/0.10	A52.0205-04
	Objective, infinity plan achromatic 10x/0.25	A52.0205-10
	Objective, infinity plan achromatic 20x/0.25	A52.0205-20
	Objective, infinity plan achromatic 40x/0.65(Spring)	A52.0205-40
	Objective, infinity plan achromatic 100x/1.25(Spring, oil)	A52.0205-100
	Objective, infinity phase contrast plan 10x/0.25 PHP	A5C.0233-10
	Objective, infinity phase contrast plan 20x/0.40 PHP (Spring)	A5C.0233-20
	Objective, infinity phase contrast plan 40x/0.65 PHP (Spring)	A5C.0233-40
	Objective, infinity phase contrast plan 100x/1.25 PHP (Spring, oil)	A5C.0233-100
	0.4x	A55.0202-01
CCD adapter	0.5x	A55.0202-04
	1x	A55.0202-02
	0.5x with dividing 0.1mm/Div	A55.0202-03
Digital Camera Adapter	For Canon EF ditigal SLR camera adapter	A55.0204-02
Photo Adapter	2.5x/4x change over photograph attachment with 10x viewing eyepiece	A55.0201-01
	4x focusing photograph attachment	A55.0201-02
	MD adapter	A55.0201-03
	PK adapter	A55.0201-04
Polarizing Attachment	Polarizer 360° rotatable	A5P.0210-1
	Analyzer	A5P.0211-1



A12.0201-A1

A12.0201-B1



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China