



CE Approval A14.0301 Trinocular Inverted Microscope 50-800x Long Working Distance

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A14.0301
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Certification: CE
- Drawtube: Binocular|Trinocular
- Illumination: Halogen Lamp
- Highlight: **inverted biological microscope,
inverted light microscope**



Product Description

50-800x Long Working Distance Trinocular Head Inverted Optical Microscope CE A14.0301

Specification	A14.0301 Biological Microscope Specification	B	T
Head	Binocular Head, inclined 30°, interpupillary distance 55-75mm	●	
	Trinocular Head, inclined 30°, interpupillary distance 55-75mm		●
Eyepiece	WF10x/18mm	●	●
	WF16x/11mm	●	●
Nosepiece	Quadruple nosepiece	●	●
Objective	Long Working Distance Plan Achromatic 4x/0.1, W.D.=17.4mm	●	●
	Long Working Distance Plan Achromatic 10x/0.25, W.D.=8.0mm	●	●
	Long Working Distance Plan Achromatic 20x/0.40, W.D.=5.0mm	●	●
	Long Working Distance Plan Achromatic 40x/0.6(S), W.D.=3.0mm	●	●
	Auxiliary objective 1.25x	●	●
	Thickness in the bottom of culature ware 1.5mm	●	●
Working Stage	Mechanical stage size 190x145mm, moving range 75x50mm, minimum vernier division 0.1mm	●	●
Focusing	Coaxial coarse & fine focusing knob, coarse focusing range 30mm, fine focusing division 0.002mm	●	●
Filter	Yellow filter Φ43mm	●	●
	Green filter Φ43mm	●	●
Light	12V 50W Halogen lamp, brightness adjustable	●	●
Optional Accessories		Item No.	
Adapter	CCD adapter	A55.0302	
	Photo adapter	A55.0301	
	Digital camera adapter	A55.0303	



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China