



Optical Binocular Stereo Zoom Microscopes WF10x A22.1219 With CE

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A22.1219
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Certification: CE
- Drawtube: Binocular
- Theory: Stereo Microscope
- Highlight: **stereo zoom microscopes,**
stereo binocular microscopes



More Images



Product Description

Optical Binocular Stereo Zoom Microscopes WF10x A22.1219 With CE

A22.1219 Stereo Microscope

Head	Binocular Head, 45° inclined
Eyepiece	WF10x/20mm
Objective	1x2x3x/1x2x4x
Working Distance	57mm
Base	With dia 70mm Glass Stage and White/Black Stage

A22.1219 Technical Data

	Magnification	View Field	Magnification	View Field	Magnification	View Field	Magnification	View Field
1X	5X	20mm	10X	20mm	15X	15mm	20X	10mm
2X	10X	10mm	20X	10mm	30X	7.5mm	40X	6.5mm
3X	15X	7mm	30X	6.7mm	45X	5.0mm	60X	4.3mm
4X	20X	5mm	40X	5mm	60X	3.5mm	80X	3.2mm



A22.1219





Eyepiece WF10/20MM





1X/2X/3X Objective



Back Side View



A22.1219



A22.1219



A22.1219



A22.1219



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China