



Industry 0.7x - 5x Monocular Zoom Video Microscope LED A21.0902

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A21.0902
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



A21.0902-C Monocular Zoom Video Microscope
With Coaxial LED Light Source, With C Mount

Product Specification

- Applications: Industry
- Certification: CE
- Drawtube: Monocular
- Illumination: LED
- Highlight: **stereo zoom microscopes,
stereo binocular microscopes**

for more products please visit us on cnoec.com

Product Description

Industry 0.7x - 5x Monocular Zoom Video Microscope LED A21.0902

A21.0902 Monocular Zoom Video Microscope

- 0.7~5x, Continus Zoom
- Objective Lens 1.0x or 1.4x, Working Distance 100mm or 105mm
- Large Black Pole Stand, With Focus Arm
- Coaxial Focus LED Light, Power 230V/12V Available

Monocular Zoom Video Microscope		A21.0902	A21.0902-C
Zoom Body	0.7~5x, Continus Zoom		
Auxiliary Objective	1.0x, Working Distance 105mm	o	
	1.4x, Working Distance 100mm		o
Stand	Large Black Pole Stand, With Focus Arm		
Adapter	Standard C-Mount		
Light Source	No Light Source	o	
	Coaxial Focus LED Light, Power 230V/12V		o



Opto-Edu (Beijing) Co., Ltd.

0086 13911110627
 sale@optoedu.com
 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China