



## Bionocular Polarizing Light Microscope For Laboratory Research CE A15.0701

### Our Product Introduction

#### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A15.0701
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



#### Product Specification

- Application: Research, Laboratory
- Certification: CE
- Drawtube: Binocular
- Illumination: Halogen Lamp
- Highlight: **student compound microscope ,  
highest magnification microscope**

## Product Description

### Bionocular Polarizing Light Microscope For Laboratory Research CE A15.0701

A15.0701 Polarizing Microscope		A	B	TR
Head	Seidentopf Binocular Viewing Head, Inclined at 30°	●	-	-
	Seidentopf Trinocular Viewing Head, Inclined at 30°	-	●	-
	Seidentopf Binocular Viewing Head, Inclined at 30°	-	-	●
Eyepiece	Eyepiece WF10×/20 with scale of cross hair	●	●	●
	Eyepiece WF10×/20 with graticule	●	●	●
	WF10×/20	●	●	●
Strain Free Objective	Infinite Plan Achromatic Objective 4x, 10x for Reflected/Transmitted system	●	●	●
	Infinite Plan Achromatic Objective 20×, 40×, 60× for Transmitted system	●	●	●
Nosepiece	Backward quadruple nosepiece, Center adjustable	●	●	●
Focusing	Coaxial coarse & fine adjustment, Range 24mm, Fine division 0.002mm	●	●	●
Stage	Revolving Round Stage Diameter Φ160mm, Graduated in 1° Increments, Minimum resolution 6' when using vernier scale	●	●	●
	Mechanical stage, Moving rang 30X40mm	●	●	●
Analyzer	Rotatable analyzer with gradation 0°-360°	●	●	●
Polarizer	360° Rotatable	●	●	●
Bertrand Lens	Bertrand Lens, focus/center adjustable sliding in/out of optical path	●	●	●
Optical Compensator	λ Slip(First class red), 1/ 4λ Slip, Quartz wedge	●	●	●
Condenser	NA0.9 Abbe condenser with Iris diaphragm & filter	●	●	●
Illumination	12V/30W Halogen lamp, Brightness adjustable(Transmitted illumination)	●	●	●
Epi-Illumination	12V/ 50W Halogen lamp with iris aperture diaphragm and iris field diaphragm, Brightness adjustable	-	-	●
Micrometer		●	●	●
Filter		●	●	●
Optional Accessories		Item No.		
Strain Free Objective	Infinite Plan Achromatic Objective 20x for Reflected system			
	Infinite Plan Achromatic Objective 40x for Reflected system			
	Infinite Plan Achromatic Objective 60x for Reflected system			
Filter	Blue			
Specimen presser				
Photo Attachment	C Mount 0.4x			
	2xSLR Camera adaptor			



A15.0701-A

A15.0701-TR



**Opto-Edu (Beijing) Co., Ltd.**

☎ 0086 13911110627

✉ [sale@optoedu.com](mailto:sale@optoedu.com)

🌐 [cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China