



Zoom Stereo Optical Microscope With Tilting Head , 1 - 4.5x , A23.1015

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE,
- Model Number: A23.1015
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Highlight: **stereo zoom microscopes,
stereo binocular microscopes**

Product Description

Zoom Stereo Optical Microscope With Tilting Head , 1 - 4.5x , A23.1015

--The microscope is among the best in its class with excellent optical performance.

--Comfortable operation with proper eyepoint

--Compact structure, long working distance, suitable for various applications.

A23.1015 Stereo Microscope		
Head	Binocular Head	
Eyepiece	EW10x	
Objective	1x-4.5x	
Zoom Ratio	1:4.5	
Work Distance	97mm	
Transmitted Illumination	6V/20W Halogen Lamp	
Incident Illumination	6V/15W Halogen Lamp with reflector	
Optional Accessories		Item No.
Head	Trinocular Head	A53.1013
Eyepiece	WF15X	A51.1021-15
	WF20X	A51.1021-20
	WF30X	A51.1021-30
Auxiliary Objective	2.0x	A52.1025-20
	0.75X	A52.1025-75
	0.63x	A52.1025-63
	0.5x	A52.1025-05
Adapter	Photo Adapter	A55.1024
	CCD Adapter 0.5x	A55.1025-05
	CCD Adapter 1.0x	A55.1025-10



www.cnoec.com



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China