



Cold Microscope Light Source 40W Microscope Accessories A56.1025

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A56.1025
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Microscope
- Certification: CE
- Product Name: Microscope Accessories
- Highlight: **microscope accessory, parts of the microscope**

Product Description

Cold Microscope Light Source 40W Microscope Accessories A56.1025


Description:

- Large Power LED Cold Light Source, 40W
- Brightness Continously Adjustable
- Wide Range Voltage 100~240V

A56.1025 LED Cold Light, 40W		
Light Source	Large Power LED 40W	
Brghtness	Continously Adjustable	
Working Life	More Than 20000 Hours	
Power	Wide Range Voltage 100~240V	
A56.1025 LED Cold Light, 40W, Optional Accessories		
	LED Cold Light Main Unit	A56.1025-1
Light Guide	Single, Rigid Flectional Tube, Length 500mm	A56.1025-2
	Dual, Rigid Flectional Tube, Length 500mm	A56.1025-3
	Ring, Flexible Tube, Length 1000mm, Size for connecting Microscope Φ60mm	A56.1025-4
Kohler	Kohler Illumination Attachment	A56.1025-5



A56.1025 LED Cold Light Source, 40W





Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China