



Binocular Head Polarized Light Microscope With Brightness Adjustable CE A15.1302

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A15.1302
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Application: Research|Laboratory
- Certification: CE
- Drawtube: Binocular
- Illumination: LED
- Theory: Polarizing Light Microscope
- Condenser: Abbe N.A.1.25
- Highlight: **student compound microscope, highest magnification microscope**



More Images



Product Description

Binocular Head Polarized Light Microscope With Brightness Adjustable CE A15.1302

Feature

A15.1302 Polarizing Microscope, with full polarizing function and competitive price, is the best economic tools for school and student use, newly developed in 2011.

- WF10x/18mm Wide Field High Eyepoint, One With Retical Micrometer
- Round Stage Dia.120mm, 360° Rotatable, Graduation 1° Increment, Minimum Resolution 6'
- Light Source 1W LED
- Polarizer: On Light Collector, 360o Rotatable
- Analyzer: Rotatable With Graduation 0~90o, Sliding Out of Optical Path
- Optical Compensator : λ Slip (First Class Red), $1/4\lambda$ Slip

Specification

	A15.1302-A	A15.1302-B
Head	Monocular Head, 30° Inclined,	Binocular Head, 30° Inclined, Interpupillary Distance 55~75mm
Eyepiece	WF10x/18mm, With Retical 0.1mm/Div	WF10x/18mm, High Eyepoint, 1 pc WF10x/18mm, With Retical 0.1mm/Div, 1 pc
Nosepiece	Backward, Quaduple Nosepiece	
Objective	Achromatic 4x,10x,40x(S),100x(O,S), 45mm DIN	
Working Stage	Round Stage Dia.120mm, 360° Rotatable, Graduation 1° Increment, Minimum Resolution 6', When Using Venire Scale Center Adjustable	
Condenser	Abbe N.A.1.25, With Dia.2~30mm Iris Diaphragm, With Filter Dia.32mm	
Focusing	Coarse & Fine Focusing Knob, Range 14mm, Minimum Division Of Fine Focusing 0.01mm	
Light	Build-in LED, 3.5V/1W, Brightness Adjustable. Polarizer: On Light Collector, 360o Rotatable Analyzer: Rotatable With Graduation 0~90o, Can Be Slided Out of Optical Path	
Polarizing	Optical Compensator : -- λ Slip (First Class Red), -- $1/4\lambda$ Slip	



A15.1302 Polarizing Microscope





Round stage, $\Phi 120\text{mm}$.360° Rotatable



Polarizer Under Condenser, 360° Rotatable



Binocular Head Interpupillary 55-75mm



Eyepiece Tube with Diopter Adjustable $\pm 5^\circ$



Tension Adjustable Focusing Knob



Round stage, $\Phi 120\text{mm}$.360° Rotatable



Condenser with Filter Holder



Rotatable analyzer with graduation 0-90°





A15.1302 Polarizing Microscope Back View



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China