



1.3M COMS 40x - 1000x Digital Optical Microscope A31.1010 For Laboratory

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE,
- Model Number: A31.1010
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Laboratory
- Certification: CE
- Drawtube: Binocular
- Illumination: Halogen Lamp
- Max Magnification: 1000X
- Theory: Digital Optical Microscope
- Highlight: **digital usb microscope,**
portable digital microscope



Product Description

1.3M COMS 40x - 1000x Digital Optical Microscope A31.1010 For Laboratory

A31.1010 Digital Laboratory Biological Microscope, Newly designed with comfortable operating with ergonomic structure, with better illumination, convenient replacement of attachment with sliding-in center adjustable condenser.

A31.1010 Digital Microscope		A	B
Head	Compensation Free Binocular Head at 30°, Interpupillary Distance 55-75mm	●	-
	Compensation Free Binocular Head at 30°, Interpupillary Distance 48-75mm	-	●
Eyepiece	High Point and Wide Field EW10/20		
Objective	Semi-plan Achromatic Objective 4x		
	Semi-plan Achromatic Objective 10x		
	Semi-plan Achromatic Objective 40x		
	Semi-plan Achromatic Objective 100x		
Nosepiece	Backward Quadruple Nosepiece		
Stage	Double Layers Mechanical Stage 140×140mm/ 75×50mm		
Condenser	N.A.1.25 Abbe Condenser		
Audio-video System	1,3M pixels CMOS digital camera with high-resolution, 1280×1024 digital camera system, USB2.0 output		
Filter	Blue Filter		
	Green Filter		
Software	Scopelmaage Advanced Micro-image Process Software		
Illumination	Halogen Lamp 6V/20W, Brightness Adjustable		
A31.1010 Digital Microscope			
Objective	Plan Achromatic Objective 4x		A52.1003-4
	Plan Achromatic Objective 10x		A52.1003-10
	Plan Achromatic Objective 20x		A52.1003-20
	Plan Achromatic Objective 40x		A52.1003-40
	Plan Achromatic Objective 100x		A52.1003-100
Illumination	Kohler Illumination		A56.1005
Dark-field Condenser	Dark-field Condenser(Dry)		A5D.1013-D
	Dark-field Condenser(Immersion)		A5D.1013-I
Polarization Set	Simple Polarization Attachment		A5P.1006
Phase Contrast Kit	Turret Phase Contrast Kit		A5C.1007
Software	Scopelmaage Analysis		A30.1004
	Scopelmaage Advanced		A30.1005





High Sensitive, Real-time Image Display by 1.3M pixel CMOS



Easy Lamp Replacement Device, just Take away the Collector



Handle for people easy to carry



Advance Software



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China