



## OPTO-EDU A15.2603-A Polarizing Microscope, Transmit Light. Binocular

Our Product Introduction

for more products please visit us on [cnoec.com](http://cnoec.com)

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A15.2603-A
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Certification: CE
- Drawtube: Binocular|Trinocular
- Illumination: Halogen Lamp|LED
- Highlight: digital trinocular microscope,  
infinity plan microscope



### More Images



Our Product

### Product Description

#### Polarizing Microscope, Transmit, Binocular

Transmit 5W LED Illumination + Center Adjustable Quadruple Nosepiece  
Non-Stress Infinity Plan Objective 4x,10x,20x,40x,60x,100x Optional  
Non-Stress Abbe N.A.1.25 Condenser, Center Adjustable With Polarizer Lens In Bottom  
Round Stage Diameter Dia.142mm. Can Be Locked And Center Adjustable

Professional Analyzer Media Unit With Bertrand Lens,  $\lambda$ ,  $\lambda/4$  Slip, Quartz Wedge.



OPTO-EDU (BEIJING) CO., LTD.

F-1501 Wanda Plaza, No.18 Shijingshan Road, Beijing 100043, China  
Tel:+8610 88696020 Fax:+8610 88696085

## A15.2603

### Polarizing Microscope, Transmit

Items	A15.2603 Transmit Polarizing Microscope Specification	A	B	Cata. No.
Head	Seidentopf Binocular Head Inclined 30°, Eyetube Rotatable 360°, Interpupillary Distance 48-75mm	●	○	A53.2621-B
	Seidentopf Trinocular Head Inclined 30°, Eyetube Rotatable 360°, Interpupillary Distance 48-75mm, , Light Split Fixed P20:E80	○	●	A53.2621-T
Eyepiece	WF10x/20mm, High Eyepoint,	●	●	A51.2602-1018
	WF10x/20mm, With Micrometer X-0.1mm*180 Div / Y-Single Line, Diopter Adjustable	●	●	A51.2605-1018D
Nosepiece	Quadruple, Backward	●	●	A54.2610-N08
Non-Stress Infinity Plan Objective	PLAN 4x/0.10, W.D.=12.1mm	●	●	A5P.2630-4
	PLAN 10x/0.25, W.D.=4.64mm	●	●	A5P.2630-10
	PLAN 20x/0.40(S), W.D.=2.41mm	●	●	A5P.2630-20
	PLAN 40x/0.66(S), W.D.=0.65mm	●	●	A5P.2630-40
	PLAN 60x/0.80(S), W.D.=0.33mm	○	○	A5P.2630-60
	PLAN 100x/1.25(S,Oil), W.D.=0.12mm	○	○	A5P.2630-100
Working Stage	Round Stage Diameter $\Phi$ 142mm, Rotatable 360°, Division 1', Vernier Division 6', Can Be Locked, Center Adjustable	●	●	A54.2602-S03
	Attachable Mechanical Stage For Polaring Microscope, Moving Range 40x40mm	○	○	A5P.2616
	Attachable Mechanical Stage For Polaring Microscope, Moving Range 30x45mm, Lecia Type, Avoid Hitting Objective Lens	○	○	A5P.2616-LK
Condenser	<b>Non-Stress</b> Abbe N.A.1.25, Iris Diaphragm, Center Adjustable With Polarizer Lens In Bottom	●	●	A56.2614-P6
Polarizing Unit	Analyzer Media Unit 0~90° Adjustable, Between Head And Body, With Built-in Betrand Lens, Removeable Out of Optical Path For Single Polarizing Observe	●	●	A5P.2621-S
	Compensator Slide, 3 Holes, With $\lambda$ Slip (First Class Red) + 1/4 $\lambda$ Slip.	●	●	A5P.2613-22a
	Quartz Wedge Slide (I-IV Class)	●	●	A5P.2615-17
	Polarizer Lens, Bulit-in Bottom Part Of Condenser, Removeable Out of Optical Path	●	●	A5P.2610-S
Transmit Illumination	5W LED Lamp Powered By Safety Low Voltage Charger, Input Wide Voltage 100V-240V, Output 5V1A, Support Power Bank Supply For Outdoor Use	●	●	A56.2651-5W
	6V/30W Halogen Lamp	○	○	A56.2613-6V30W
Filter	Dia.45mm, Blue (Built in)	○	○	A56.2616-45B
	Dia.45mm, Green	○	○	A56.2616-45G
	Dia.45mm, Yellow	○	○	A56.2616-45A
	Dia.45mm, Neutral	○	○	A56.2616-45N
Adapter	C-Mount 0.5x, Focus Adjustable	○	○	A55.2601-05
	C-Mount 0.75x, Focus Adjustable	○	○	A55.2601-75
	C-Mount1.0x, Focus Adjustable	○	○	A55.2601-10
	Digital Eyepiece Camera Adapter, Dia. 23.2mm/25mm	○	○	A55.2610
Package	Carton Size 570×275×450mm, 1 Unit/Ctn			
Weight	Gross Weight: 8.6 kgs, Net Weight: 7.5 kgs			

Note: "●"In Table Is Standard outfits, "○" Is Optional Accessories.



OPTO-EDU (BEIJING) CO., LTD.

F-1501 Wanda Plaza, No.18 Shijingshan Road, Beijing 100043, China  
Tel:+8610 88696020 Fax:+8610 88696085

## A15.2603

### Polarizing Microscope, Transmit



OPTO-EDU (BEIJING) CO., LTD.

F-1501 Wanda Plaza, No.18 Shijingshan Road, Beijing 100043, China  
Tel:+8610 88696020 Fax:+8610 88696085

**A13.2003**

## **Polarizing Microscope, Transmit**



**Condenser & Light Source**  
Abbe Condenser N.A.1.25, Iris Diaphragm, Rack & Pinion Up/Down, With Filter Holder  
3W-LED Illumination Systems



**Coaxial Coarse & Fine Focusing**



**Wire Winding Device**




**Holding Hole On Back**

**Powered By Safety Low Voltage Charger**  
Input Wide Voltage 100V-240V, Output 5V1A, Support Power Bank Supply For Outdoor Use



## Product Size(mm)



## Product Multi-angle



## Product Accessories



**A55.2601-05**  
C-Mount 0.5x  
Focus  
Adjustable



**A55.2601-75**  
C-Mount 0.75x  
Focus  
Adjustable



**A55.2601-10**  
C-Mount 1.0x  
Focus  
Adjustable




**A55.2610**  
Digital Eyepiece  
Camera Adapter  
Dia. 23.2mm/25mm



**A5P.2616-LK** Attachable Mechanical  
Stage, Lecia Type



**A5P.2616** Attachable Polarizing  
Mechanical Stage, Moving Range 40x40mm



**Opto-Edu (Beijing) Co., Ltd.**



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China