



Lab Student Achromatic Binocular Compound Microscope With LED Light

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1137
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Highlight: binocular biological microscope, biological compound microscope



More Images



for more products please visit us on cnoec.com

Product Description

Lab Student Achromatic Binocular Compound Microscope With LED Light

A11.1137 Student Microscope,40x-400x

- 40x~1000x Magnification
- LED Light Source
- Full Metal Body for Steday Observation
- Coaxial Coarse & Fine Focusing Knob

A11.1137 Student Biological Microscope

Head	Binocular Head, Diopter Adjustable
Eyepiece	WF10x/18mm
Objective	Achromatic 4x,10x,40x,100x (Spring Oil)
Nosepiece	Backward Quaduple Revolving Nosepiece
Stage	Double Layers Mechanical Working Stage
Condenser	Abbe Condenser N.A.1.25, With Iris Diaphargm
Focusing	Coaxial Coarse & Fine Focusing Knob
Illumination	LED Light Source, Brightness Adjustable




A11.1137 Biological Microscope





A11.1137 Biological Microscope





A11.1137 Biological Microscope



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China