



Infinity Industry Monocular / Trinocular DIC Metallurgical Optical Microscope A13.0217

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A13.0217
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month





Product Specification

- Highlight: **reflected light microscope, halogen lamp microscope**



for more products please visit us on cnoec.com

More Images



Product Description

Infinity Industry Monocular / Trinocular DIC Metallurgical Optical Microscope A13.0217

A13.0217 DIC Metallurgical Microscope, used for biology medicine agriculture and industry area widely. It is the ideal instrument in medical treatment, teaching demonstration and scientific research.

- High Quality Infinity Optical System
- Monocular, Trinocular Head Optional
- DIC Slide For 5x/10x/20x Objective
- Strainfree Infinity LWD Plan Achromatic 5x,10x,20x Optional
- Polarizing System, Switch In/Out Of Light Path

A13.0217 DIC Metallurgical Microscope		A	B
System	Infinity Optical System, Total Magnification 100x Mechanical Tube Length ∞ Objective Conjugate Distance∞		
Head	Monocular Trinocular Head, 30° Inclined, Light Separate Ratio 80:20	o	o
Eyepiece	WF10x/18mm, Eyepiece Tube Dia.23.2mm, Parfocal Distance 10mm		
Objective	Strainfree LWD Infinity Plan Achromatic PL10x/0.25, W.D.20.2mm		
Working Stage	Double Layer Mechanical Stage 180*145mm, Moving Range 35*30mm		
Focusing	Coaxial Coarse & Fine Focus Knob, Tension Adjustable, Minimum Division 0.002mm, Height Adjustabout 133mm, Max Sample Height 130mm		
DIC	DIC Slide Match 5x/10x/20x Objectives, Horizontal Adjustable		
Polarizing	Polarizer, 360° Rotatable, Switch In/Out Analyser, Switch In/Out		
Light Source	Coaxial Illumination, High Light LED, Brightness Adjustable		
Power	Wide Range Voltage 85~265V, 50/60Hz		
Optional Accessories			
Eyepiece	WF15x/11mm		A51.0203-16A
	WF20x/11mm		A51.0203-20A
Objective	Strainfree LWD Infinity Plan Achromatic PL5x/0.1,W.D.18.2mm		A5M.0235-5A
	Strainfree LWD Infinity Plan Achromatic PL20x/0.4,W.D.8.8mm		A5M.0235-20
Adapter	CCD Adapter 0.4x		A55.0202-01
	CCD Adapter 0.5x		A55.0202-04
	CCD Adapter 1.0x		A55.0202-02
	CCD Adapter 0.5x With Dividing 0.1mm/Divv		A55.0202-03
	Digital SLR Camera Adapter for Canon EF Series		A55.0204-02



A13.0217-B Trinocular Metallurgical Microscope, DIC

A13.0217-B Monocular Metallurgical Microscope, DIC

www.cnoec.com




A13.0217-A Use With Digital SLR Camera



High Quality DIC Objective, DIC Slide



Epi-Illumination System, Coaxial Illumination Using High Brightness White LED, Brightness Adjustable, Power Box 85~265V




Metal Base 300*250mm, Mechanical Stage 160*140mm, Moving Range 35*30mm, X/Y Coaxial Moving Knob On Stage




A13.0217-A Use With A59.0204 Digital Camera



A13.0217-A Use With A59.2208 Digital Camera



A13.0217 Metallurgical Microscope System Diaphragm



A13.0217 Metallurgical Microscope System Size



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China