



Video 40x Handheld Portable Digital Microscope 18mm Eyepiece

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A31.0921
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month




Product Specification

- Function: Handheld Video Microscope
- Head: Monocular Head, 30° Inclined
- Eyepiece: PL10X/18mm
- Nosepiece: Triple Nosepiece
- Objective: Achromatic 4x/0.10, 10x/0.25, 40xs/0.65
- Condenser: NA0.65 Condenser And Iris Diaphragm
- Magnification: 40x~400x
- Light Source: LED Illumination (3 Bulbs)
- Highlight: **opto edu handheld portable digital microscope, 40x handheld portable digital microscope, 40x handheld digital microscope**



More Images



Product Description

WIFI Digital Biological Microscope, 2.0M

Monocular Head Built-in WIFI Camera 2.0MP

Support Android Tablet

Achromatic 4x/0.10, 10x/0.25, 40xs/0.65

N.A.0.65 Condenser With Iris Diaphragm

LED Illumination 3 Bulbs Brightness Adjustable



OPTO-EDU (BEIJING) CO., LTD.

OPTO-EDU

F-1501 Wanda Plaza, No.18 Shijingshan Road, Beijing 100043, China
Tel:+8610 88696020 Fax:+8610 88696085

A31.0921-2M WIFI Digital Microscope



Product Details



A31.0921	
Magnification	40x~400x
Head	Monocular Head, 30° Inclined Built-in WIFI camera, 2MP. With Pad holder
Eyepiece	PL10X/18mm
Nosepiece	Triple Nosepiece
Objective	Achromatic 4x/0.10, 10x/0.25, 40xs/0.65
Stage	Metal Plain Stage 110*120mm, With Slide Clip
Focusing	Separate Coarse & Fine Focusing. Coarse Focusing Range 8.7mm, Fine Focusing Range 1.2mm
Condenser	NA0.65 condenser and iris diaphragm
Light Source	LED Illumination (3 bulbs)



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China