




Stereo Usb Portable IPS Screen Hand Held Digital Microscope Dual Lens

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A36.5102
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Function: Wireless Digital Microscope Handheld
- LCD Screen: "7"" IPS Screen, With Touch Button To Take Photo & Video, One Key Switch Between Main Body Camera & Portable USB Camera"
- Handheld Camera: "1.3M Digital Camra, Zoom Lens 1.5x~25x, Output To 7"" LCD Screen 8 LEDs Light, Brightness Adjustable"
- Main Body Camera: "2.0M Digital Camra, Zoom Lens 5x~35x, TF Card, USB Output Coaxial LED Ring Light, Brightness Adjustable"
- TF Card Output: Support Max 32G, To Take Photo & Video
- USB Output: To Win/Mac, With Professional StrangeView Measuring Software
- Stand: Pole Stand, With Focusing Knob, Size 180x170mm
- Size: 180*170*322mm
- Power Supply: DC 5V2A Adapter, Can Be Powered By USB Power Bank For Outdoor Use



More Images



Trinocular Head With WF10x/20mm Eyepiece

Product Description

7" LCD Stereo + USB Portable Dual Lens Digital Microscope, 2.0M+1.3M


Brand New Design Dual Lens LCD + Portable Digital Stereo Microscope For Different Purposes

0.7x~4.5x Zoom Stereo Trinocular Microscope With WF10x/20mm Eyepiece

7.0" LCD, Main Body 2.0M Camera, 5x-35x Zoom For Real Time Clear Watching

Portable USB2.0 1.3M Digital Microscope, 1.5x-25x Zoom For Large Size Objective Watch on LCD

USB Output To Win/Mac, With Professional StrangeView Measuring Software



7" IPS Screen

2.0M Digital Camera, Support Max 32G TF Card, USB Output To Win/Mac



Head

Trinocular

Eyepiece

WF10x/20mm



Excellent Optical Lens

Ultra-clear Focusing System Increase The Image Quality & Supply Wide Field Of View





LED Ring Light

Main Body Illumination Coaxial LED Ring Light,
Brightness Adjustable



A36.5102 7" LCD Stereo + USB Portable Dual Lens Digital Microscope, 2.0M+1.3M	
LCD Screen	7" IPS Screen, With Touch Button To Take Photo & Video, One Key Switch Between Main Body Camera & Portable USB Camera
Handheld Camera	1.3M Digital Camra, Zoom Lens 1.5x~25x, Output to 7" LCD Screen 8 LEDs Light, Brightness Adjustable
Main Body Camera	2.0M Digital Camra, Zoom Lens 5x~35x, TF Card, USB Output Coaxial LED Ring Light, Brightness Adjustable
TF Card Output	Support Max 32G, To Take Photo & Video
USB Output	To Win/Mac, With Professional StrangeView Measuring Software
Stand	Pole Stand, With Focusing Knob, Size 180x170mm
Size	180*170*322mm
Power Supply	DC 5V2A Adapter, Can Be Powered By USB Power Bank For Outdoor Use
Package	Color Box + Foam, 1 pc/1 Box, With Software Disc, DC Adapter, Instruction Manual



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China