



A11.0205 Biological Microscope Binocular Microscopes WF 10X

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.0205
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Highlight: binocular biological microscope,
biological compound microscope

Product Description

Lab Achromatic Biological Microscope Binocular Microscopes A11.0205

This microscope is the conveniently operated optical instruments with excellent functions and new pattern it is very suitable for usage in medical test, teaching and microobservation in hospitals and schools.

--A, C, with coaxial coarse/fine focusing system, with up stop.

--B, D, with coaxial coarse/fine focusing system, with tensional adjustable and up stop

--A1,B1,C1,D1 is monocular head, A2,B2,C2,D2 is binocular head

--12V 20W halogen lamp with brightness control

pecification	A11.0205-A	A11.0205-B	A11.0205-C	A11.0205-D
Eyepiece	WF 10X(Φ18mm)			
Objective	Achromatic 4X/0.10			
	Achromatic 10X/0.25			
	Achromatic 40X/0.65(Spring)			
	Achromatic 100X/1.25(Spring, oil)			
Head	Compensation binocular head inclination 30°			
Focus System	Coaxial coarse/fine system focus, Minimum division of fine focusing: 4um			
	With Up stop	With tensional adjustable and up stop	With Up stop	With tensional adjustable and up stop
Nosepiece	Quadruple Frontward ball bearing outer locating		Quadruple Backward ball bearing inner locating	
Stage	Double layer mechanical, size 125x115mm, moving range 70x30mm			
Condenser	N.A.1.25 Rack & pinion adjustment			
Filter	Blue filter			
	Frosted filter			
Collector	For halogen lamp			
Light source	12V 20W halogen lamp with brightness control			



Optional Accessories		Item No.
Eyepiece	Wide field WF16x/11mm	A51.0203-16A
	Dividing 10x/17mm	A51.0205-10
Head	Monocular, inclination of 45°	A53.0201-3
	Sliding binocular, inclination of 30°	A53.0203-121
	Dual viewing, inclination of 30°	A53.0202-1
Nosepiece	Triple nosepiece(Frontward ball bearing outer locating)	A54.0201
Stage	Fixed with slide moving holder, size 120x120mm, move 70x30mm	A54.0205-05
	Fixed, size 120x120mm	A54.0205-01
Filter	Green filter	A56.0209-2
	Yellow filter	A56.0209-3
Light Source	Fluorescent lamp 5W, 110V/220V	A56.0202
	Incandescent lamp, 110V/220V	A56.0203
	LED Illumination, high brightness white	A56.0204
	Transformer for LED, Input 230V, 50/60Hz, Output 9V500mA	A56.0205-1
	Transformer for LED, Input 110V, 50/60Hz, Output 9V500mA	A56.0205-2
	Reflecting mirror, plane-concave reflecting mirror with stand	A56.0206
Polarizing Set	Polarizer, 360° rotatable	A5P.02010-1
	Analyzer	A5P.02011-1



A11.0205-A2



A11.0205-C2



 0086 13911110627  sale@optoedu.com  cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China