



## Binocular Metallurgical Halogen Lamp Microscope , Plan Achoromatic Optical System

Our Product Introduction

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A13.0902
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Highlight: reflected light microscope ,  
bright field microscope



for more products please visit us on [cnoec.com](http://cnoec.com)

## Product Description

### Binocular Metallurgical Halogen Lamp Microscope , Plan Achoromatic Optical System

Specific ation	A13.0902	A B
Optic Syste m	Plan Achoromatic Optical System	• •
Head	Compensation free Binocular,30°inclined Compensation free Trinocular,30°inclined	• •
Eyepe ce	PL 10X High Eye-point plan eyepiece	
Object ives	195 Plan achromatic metallurgical 4X,10X,40X	
Nosepi ece	Triple or quadruple nosepiece, setting accuracy 0.03mm	
Stage	Rectangular and movable mechanical stage with low coaxial adjusting(right and left hand for option),140*132mm(accompanying 180*145mm),working distance:75*50mm	
Focus	Coaxial focus system with tension adjustment,with up-limited coarse focus, Stage range 25mm, precision 0.002mm	
Illumin ation Unit	Adopted 90-240V voltage,Kohler illuminator system, the 6V/30W Halogen bulb,adjustable filament center and adjustable brightness,built-in field diaphragm.	
Optional Polariz er	Accessories Adjusting range:plusminus $\pm 25^\circ$	



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



[sale@optoedu.com](mailto:sale@optoedu.com)



[cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China