



## Trionocular Polarizing Light Microscope For Laboratory Research CE

### Our Product Introduction

#### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A15.0702
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



#### Product Specification

- Highlight: **student compound microscope ,  
highest magnification microscope**

## Product Description

### Trinocular Polarizing Light Microscope For Laboratory Research CE

A15.0702 Polarizing Microscope		A	B
Head	Sliding Trinocular Viewing Head, Inclined at 30°	●	●
Eyepiece	Eyepiece WF10×/18 with scale of cross hair	●	●
	Eyepiece WF10×/18 with graticule	●	●
	WF10×/18	●	●
Objective	Strain Free Achromatic Objective 4×,10×	●	●
	Strain Free Achromatic Objective 25×, 40×, 63× for Transmitted system	●	●
	Strain Free Achromatic Objective 25×, 40×, 63× for Reflected system	●	●
Nosepiece	Backward quadruple nosepiece, Center adjustable	●	●
Focusing	Coaxial coarse & fine adjustment, Range 28mm, Fine division 0.002mm	●	●
Stage	Revolving Round Stage Diameter Φ160mm, Graduated in 1° Increments, Minimum resolution 6' when using vernier scale	●	●
	Mechanical stage, Moving rang 30X40mm	●	●
Analyzer	Rotatable analyzer with gradation 0°-90°	●	●
Polarizer	360° Rotatable	●	●
Bertrand Lens	Bertrand Lens, sliding in/out of optical path	●	●
Optical Compensator	λ Slip (First class red), 1/ 4λ Slip, Quartz wedge	●	●
Condenser	NA0.1.25 Abbe condenser with Iris diaphragm & filter	●	●
Illumination	12V/20W Halogen lamp, Brightness adjustable (Transmitted illumination)	●	●
	12V/ 50W Halogen lamp Brightness adjustable (Reflected illumination)	-	●
Micrometer		●	●
Optional Accessories		Item No.	
Objective	Lamp house adjustment objective		
Filter			
Specimen presser			
Photo Attachment	C Mount 0.5xC		
	C Mount 1xC		



**Opto-Edu (Beijing) Co., Ltd.**

☎ 0086 13911110627

✉ [sale@optoedu.com](mailto:sale@optoedu.com)

🌐 [cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China