



40X - 1000X Trinocular Compound Microscopes Fluorescence A16.0301

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A16.0301
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Highlight: epi-fluorescence microscope, epifluorescence microscope

Product Description

40X - 1000X Trinocular Compound Microscopes Fluorescence A16.0301

Specification	A16.0301 Biological Microscope Specification	B	T
Head	Binocular Head, inclined 30°, interpupillary distance 55-75mm	•	
	Trinocular Head, inclined 30°, interpupillary distance 55-75mm		•
Eyepiece	WF10x/18mm	•	•
Nosepiece	Quadruple nosepiece	•	•
Objective	Achromatic 4x/0.1	•	•
	Achromatic 10x/0.25	•	•
	Achromatic 40x/0.65(S)	•	•
	Achromatic 100x/1.25(O,S)	•	•
Working Stage	Mechanical stage size 140x140mm, moving range 70x50mm, minimum vernier division 0.1mm	•	•
Condenser	N.A.1.25 Abbe condenser, up/down range 15mm	•	•
Focusing	Coaxial coarse & fine focusing knob, coarse focusing range 20mm, fine focusing division 0.002mm	•	•
Fluorescent Unit	A57.0302 Epi-Fluorescent Unit --B(Blue), G(Green), U(Ultraviolet), V(Violet) exciting light filter system --100W super high pressure spherical mercury lamp --Ultraviolet light resist screen	•	•
Filter	Blue filter	•	•
	Yellow filter	•	•
	Ground glass	•	•
Light	6V 20W Halogen lamp, brightness adjustable	•	•
Optional Accessories		Item No.	
Eyepiece	WF16x	A51.0301-16	
	WF20x	A51.0301-20	
Objective	Achromatic 20x (S)	A52.0301-20	
	Achromatic 60x (S)	A52.0301-60	
	Plan 4x	A52.0303-04	
	Plan 10x	A52.0303-10	
	Plan 20x(S)	A52.0303-20	
	Plan 40x(S)	A52.0303-40	
	Plan 60x(S)	A52.0303-60	
	Plan 100x(O,S)	A52.0303-100	
Adapter	CCD adapter	A55.0302	
	Photo adapter	A55.0301	
	Digital camera adapter	A55.0303	
Condenser	Dark field condenser, dry	A57.0310-01	
	Dark field condenser, immersion	A57.0310-02	
Fluorescent	A57.0371 Epi-Fluorescent Unit --B(Blue), G(Green) exciting light filter system --100W super high pressure spherical mercury lamp --Ultraviolet light resist screen	A57.0371-01	



A16.0301-B

A16.0301-T



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China