



Industrial Stereo Optical Microscope With Double Arm Pole Stand

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE,
- Model Number: A23.0906-S2
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification


- Highlight: stereo zoom microscopes,
stereo binocular microscopes

Product Description

Industrial Stereo Optical Microscope With Double Arm Pole Stand

- 0.7~4.5x Zoom Lens, with WF10x/20mm Eyepiece,
- Total Optical Magnification 7~45x
- Double Arm Universal Pole Stand

A23.0906-S2 Zoom Stereo Microscope, 0.7x-4.5x Specification		A23.0906-S2	A23.0906-S2T
Head	Binocular, 45° Inclined, 360° Rotatable, +/- 5° Diopter Adjustable, Interpupillary Adjustment 54~75mm	o	
	Trinocular, 45° Inclined, 360° Rotatable, +/- 5° Diopter Adjustable, Interpupillary Adjustment 54~75mm		o
Zoom Range	0.7~4.5x		
Zoom Ratio	1:6.5		
Eyepiece	WF10x/20mm, Both with Eye-cup		
Magnification	7~45x, Up to 3.5~180x with Auxiliary Lens		
Working Distance	100mm, With Standard 1x Lens		
Field of View	1~95.2mm		
Stand	Double Arm Universal Pole Stand STL2, Height 380mm, Focusing Pole Height 140mm, Level Move Range 235mm, Base 210*260*20mm,		
Illumination	No Light Source		
A23.0906-S2 Zoom Stereo Microscope, 0.7x-4.5x Optional Accessories			
Eyepiece	WF15x/15mm	A51.0922-15	
	WF20x/10mm	A51.0922-20	
Auxiliary Lens	0.5x	A52.0927-05	
	0.75x	A52.0927-75	
	1.5x	A52.0927-15	
	2.0x	A52.0927-20	



A23.0906-S2 Zoom Stereo Microscope 0.7~4.5x, With Universal Stand



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China