



Nail Checking Fluorescence Digital Optical Microscope Video Biological

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE,
- Model Number: A31.0210
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Highlight: **digital usb microscope, optical microscope with digital camera**



More Images



for more products please visit us on cnoec.com

Product Description

Nail Checking Fluorescence Digital Optical Microscope Video Biological

A31.0210 Nail Fold Microcirculation Checking Microscope, also named as " **Microcirculate Checkup Microscope**" is a new type of optical & electronic instrument, equipped with achromatic objective, new type LED illumination, high definition color sensor. It is mainly applied to nail fold microcirculation checking, no any hurt and side-effect in the checking process. It is the ideal instrument for disease forepart diagnosis, health care etc.

	A31.0210-A	A31.0210-B
Objective	Achromatic 4x/0.10	
Color Sensor	Video Output	USB Output
	1/3"CMOS, 380 TV lines Automatic white balance	1/3"CMOS, 1.3M USB2.0
Power Adapter	Input 220V/110V, output 9V	
Movable Stage	Mechanism moving stage, size 130x118mm, moving range 64x28mm	
Video Cable	Standard video cable, length 1.8m	
LED Illumination	High quality LED, brightness and illuminate's angle adjustable	
Optional Accessories		Item No.
Objective	Plan achromatic 5X/0.12	A52.0203-5

www.cnoec.com



A31.0210 Nail Fold Microcirculation Checking Microscope



www.cnoec.com



A31.0210 Nail Fold Microcirculation Checking Microscope, Back View



Including Video Output Cable, Power Supply



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China