



Metallurgical Optical Microscope Infinity System Sliding Trinocular A13.1304

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A13.1304
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Highlight: reflected light microscope ,
halogen lamp microscope

Product Description

Infinity System Sliding Trinocular Metallurgical Optical Microscope A13.1304

A13.1304 Metallurgical Microscope, with high quality component design, simple and flexible combination and strong metallographic detecting function, it is the best choice for large size sample observation

- Compact and Flexible Structure, Adapting Various Surroundings
- Get Clear and Sharp Image Under Various Kind of Observations
- Optical Components Are All Specially Handled and Coated With Special Film
- New Design of Long Working Distance Professional Metallurgical Objectives

A13.1304 Metallurgical Microscope	
Head	Sliding Trinocular Head, 45° Inclined, 360° rotatable, interpupillary distance: 55-75mm;
Eyepiece	Wide Field Plane-Scope: WF10X/18mm
Nosepiece	Quadruple Nosepiece
Infinity Objective	45mm DIN Achromatic 10x, 100x (O.S)
	Plan objective of infinite correction system: PC10X, PC20X, PC40X, PC80X;
Working Stage	Double Layer Mechanical Stage, Size 190*150mm, Moving Range X-Y 90*60mm
Condenser	Abbe N.A.1.25, With Dia.2-Dia30mm Iris Diaphragm,
Focusing	Coaxial coarse and fine focusing adjustable mechanism, step-by-step variable gear, rack and pinion mechanism, with mechanism for adjusting focusing torque and controlling the upper limit of stage lift. Focusing Range: 36mm, precision: 0.002mm/scale;
Polarizer	With polarizer and analyzer. The polarizer can sliding out of optical path;
Filters	Five stop and aperture stop, Blue, green filters;
Light Source	Vertical illumination system for reflected illumination, 12V/50W halogen lamp, light adjustable continuously, bottom light system with extra light LED for transmitted light, light adjustable continuously;
Gross Weight	12KG
Package	1set/carton
Dimension	38*20*46cm





0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China