



Stereo Optical Microscope With Tilting Binocular Head 1x 3x , A22.1307

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE,
- Model Number: A22.1307
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Product Name: Stereo Microscope
- Head: Binocular
- Stand: Pole Stand
- Focus: Rack And Pinion Coarse Focusing
- Illumination: Transmit And Reflect LED, Brightness Adjustable
- Working Distance: 100mm
- Highlight: **stereo zoom microscopes,**
stereo binocular microscopes



More Images



Product Description

Stereo Optical Microscope With Tilting Binocular Head 1x 3x , A22.1307

Introduction

A22.1307 is the basic style stereo microscope whose price is affordable. Stereo microscope is also known as dissecting microscope, it is widely used in biology, medicine, agriculture, forestry, industry and marine life etc. A22.1307 has top and bottom light source, both can adjust brightness, and pole stand with spring clips is convenient for users.

Specification

A22.1307 Stereo Microscope		13	12	24
Head	Binocular Head, Inclined 45°, Interpupillary Distance 55-75mm, With A Diopter Adaptor On The Left Eyepiece Tube			
Eyepiece	WF10x/18mm			
Objective	1x/3x	√		
	1x/2x		√	
	2x/4x			√
Working Distance	100mm			
Stage	Ø95mm, 1 PC Frosted Glass Stage (For Bottom Illumination), 1 PC Plastic Black/White Plate (For Top Illumination)			
Focus	Rack And Pinion Coarse Focusing, Focusing Range 40mm			
Stand	Pole Stand, Base Size 21x14x5cm, With Spring Clips			
Illumination	Transmit & Reflect LED, Both Brightness Adjustable			
Power Supply	110~240V Wide Range Power Supply			

Product Details



01 BINOCULAR HEAD

45° inclined, with a diopter adaptor on the left eyepiece tube & WF10x/18mm eyepiece.

02 OBJECTIVE

Three different objectives for choosing (1x/3x, 1x/2x, 2x/4x) according to your needs.



03 POLE STAND

Base size 21x14x5cm, with spring clips.

04 POWER SUPPLY

110~240V, internal power supply with wide range.



05 LIGHT SOURCE

Top LED light with two adjustment Knobs, 1 piece plastic black/white plate for top illumination.


Bottom LED Light, 1 piece frosted glass stage for bottom illumination.

Stage Ø95mm.



Product Application

Stereo microscope, also known as "solid microscope" or "operation and dissecting microscope", is a kind of stereoscopic microscope. Stereo microscope is widely used in materials such as macroscopic surface observation, failure analysis, fracture analysis and other industries field. As a kind of visual equipment with three-dimensional, it plays an important role in biology, medicine, agriculture, forestry, industrial and marine biological industry.



Packaging & Shipping

Package	4 Sets/Carton, 67x40x41 CM
Gross Weight	12KG
Packing	Polyfoam Protection And Strong Carton.
Shipping Time	5~20 Days After Receive The Deposit.



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China