



0.7x - 4.5x Track Stand , No Light Source Zoom Stereo Microscope A23.1302

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE,
- Model Number: A23.1302
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Highlight: stereo zoom microscopes,
stereo dissecting microscope



More Images



Product Description

0.7x - 4.5x Track Stand , LED Light Source Zoom Stereo Microscope A23.1302

- 0.7~4.5x Zoom Lens, with WF10x/22mm Eyepiece,
- Total Optical Magnification 7~45x
- Large base with Track Stand

A23.1302 Zoom Stereo Microscope		
Head	Binocular Head, Interpupillary Distance:55-75mm,45o Inclined,360o Rotatable ,with a diopter adaptor on the both eyepiece tube	
Eyepiece	W10x/22mm	
Objective	0.7x-4.5x	
Work Distance	100mm	
Light Source	Pole Stand without light source	
Focus	Rack and Pinion Focusing,with hand wheels for adjusting focusing torque. Focusing Range:75mm;	
Optional Accessories		Item No.
Auxiliary Objective	2.0x	
	0.75X	
	0.5x	
	1.5X	



A23.1302 Stereo Microscope





A23.1302 Eyepiece WF10X/22mm



A23.1302 Focusing



A23.1302 Tension Knob



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China