



Trinocular Zoom Stereo Optical Microscope 0.7x - 4.5x LED Light Source , A23.1303

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE,
- Model Number: A23.1303
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Highlight: stereo dissecting microscope, stereo binocular microscopes



More Images



Product Description

Trinocular Zoom Stereo Optical Microscope 0.7x - 4.5x LED Light Source , A23.1303

- 0.7~4.5x Zoom Lens, with WF10x/22mm Eyepiece,
- Total Optical Magnification 7~45x
- Large base with Track Stand

A23.1303 Zoom Stereo Microscope	
Head	Binocular Head, Interpupillary Distance:55-75mm,45o Inclined,360o Rotatable ,with a diopter adaptor on the both eyepiece tube
Eyepiece	W10x/22mm
Objective	0.7x-4.5x
Work Distance	100mm
Light Source	Track Stand with Transmitted Illumination and Incident Illumination 3W LED,brightness adjustable
Base	290X245X60mm
Focus	Rack and Pinion Focusing,with hand wheels for adjusting focusing torque. Focusing Range:75mm;
Package	1 Set/Carton
Dimension	37X32X57cm
Net Weigh	5kgs
Optional Accessories	
	2.0x
Item No.	

Auxiliary Objective	0.75X	
	0.5x	
	1.5X	
Light	External Ring Light	



A23.1303 Zoom Stereo Microscope



A23.1303 Zoom Stereo Microscope





A23.1303 Zoom Stereo Microscope



A23.1303 Trinocular Head



A23.1303 Eyepiece WF10X/22mm





A23.1303 Focusing



A23.1303 Tension Knob



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627 ✉ sale@optoedu.com 🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China