



## WIFI Digital Microscope A34.5002

### Our Product Introduction

#### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A34.5002
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



#### Product Specification

- Highlight: Portable Digital Microscope, hand held microscope



#### More Images



## Product Description

### WIFI Microscope 200X 1.0M Digital Microscope iPad/Android/Win/Mac A34.5002

A34.5002 WIFI Microscope, iPhone/iPad/Pc/Android, 10X-200X

- Magnification 10x~200x
- WIFI Transmission Rate 150Mbps
- Light Source 6 LEDs, Brightness Adjustable
- Support iPhone/iPad/Pc/Android/Mac Book
- It can work on Windows/Mac/iOS (iPhone, iPad)/Android (smart phone, Tablets)
- It has measurement function (by software) when working on Windows/Mac
- It comes with polarizer (which can reduce reflection when used to observe shiny objects like jewelry, metal etc)
- Professional Stand, very suitable for observation

A34.5002 WIFI Microscope, iPhone/iPad/Pc/Android, 10X-200X	
Resolution	1.0M, True Resolution
Senser	1/4" CMOS, Color
Magnification	10x~200x
Video Resolution	1280x720 (720P), 640x480 (VGA)
Video Format	H.264
Frame Rate	Wireless 720P: 15fps, VGA: 25fps
Photo Resolution	1280x720 (720P), 640x480 (VGA)
Photo Format	JPG
Focus Range	0~40mm, Manual Focus
Snapshot	By Software
PC Interface	USB2.0
IOS/Android Interface	WIFI Direct
Transmission Range	20M, Open Sight
Transmission Rate	150Mbps
Power Supply	DC5V, 1A, by Mini USB
Battery	Rechargeable Battery, Lithium 1300mAh, Can Be Charged By Adaptor/USB
Battery Working Time	Around 75 Minutes
Operation Temperature	0~65°
Light Source	6 LEDs, Brightness Adjustable
Power Source	Rechargeable Li-ion Battery Working Time: 90mins, Charging Time: 100mins Charging Power: 5V DC
Operation System	IOS (5.1 or Later) Android Windows XP/Win7/Win8 Mac Book (10.7 or Later)
App	IOS: WifiCapture, Search And Download From App Store Android: WifiCapture, Install from CD, or Search Download From <a href="https://play.google.com/store">https://play.google.com/store</a>
OSD Language	English, Deutsch, Spanish, Russian, Italian, French, Portuguese, Dutch, Polish, Japanese, Korean, Chinese
Bundle Software For PC	PortableCaptureHD (with Measurement)
Measureement	Yes, For Windows and Mac Book
Accessories	Professional Stand, Adapter, USB Cable, Softare CD, Calibration Charge
Dimension	Microscope: 154mm (L) x 44mm (R) Professional Stand: 165mm(L) x118mm(W)x12mm (H)
Weight	130g
Packaging	Gift Box, 6pcs/carton 9.0kgs/carton, 44x42x35cm





**A34.5002 WIFI Microscope, iPhone/iPad/Pc/Android, 10X-200X**



[www.cnoec.com](http://www.cnoec.com)



[www.cnoec.com](http://www.cnoec.com)





A34.5002 WiFi Microscope, With Stand Height Adjustable



A34.5002 WiFi Microscope, Packing List



**Opto-Edu (Beijing) Co., Ltd.**

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China