



Desk Top Polariscope 360 Degree Platform Rotatable Jewelry Microscope A24.6331 - A

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE,
- Model Number: A24.6331-A
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Highlight: laboratory binocular microscope,
gemological microscope



More Images



Product Description

Desk Top Polariscope 360 Degree Platform Rotatable Jewelry Microscope A24.6331 - A

- 1, 49mm Calibre
- 2, Warm White LED Light
- 3, Voltage:DC3V(110V-240V)

A24.6331-A Desk Top Polariscope	
Polarizing Lens	49mm Big Caliber
Sample Platform	360° Rotatable
Illumination	Warm White LED
	Imitation Incandescent Light
Power Supply	Attachable DC3V (110V-220V)



A24.6331-A
Desk Top Polariscope



49mm Big Caliber Polarizing Lens
Convenient For The Interference
Pattern Observing



360° Rotatable Sample Platform
Warm White LED Illumination
Imitation Incandescent Light



Glass Ball Under Polarizing Observation



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China