



Calcite Dichroscope Observing for Pleochromatism Of Gems A24.6351

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE,
- Model Number: A24.6351
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Highlight: laboratory binocular microscope,
gemological microscope



Product Description

Calcite Dichroscope Observing for Pleochromatism Of Gems A24.6351

It Is Suitable For Observing The Pleochromatism Of Gems

A24.6351 Calcite Dichroscope	
A24.6351-A	A24.6351-B
Handheld Type	Desktop Type (With Stage)
Excellent Iceland Spar	
Full Black Main Body Avoid Reflecting Interference, Length 53mm, Diameter 16.5-18mm	
Big Eyepiece	
Standard Rectangle Window	
Without Stage	With Stage, Gemstone Can Be Adjusted, 360° Rotatable, Up & Down, Backforward & Frontward Adjustable
Leather Bag Package	



Leather Bag

Objective

Eyepiece

A24.6351-A
Handheld Calcite Dichroscope



Rotatable Knob

A24.6351-B
Calcite Dichroscope with Stage



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China