



## 100 - 400x Digital Portable Metallurgical Optical Microscope A31.2501

Our Product Introduction

for more products please visit us on [cnoec.com](http://cnoec.com)

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A31.2501
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Highlight: reflected light microscope,  
halogen lamp microscope



### More Images



### Product Description

#### 100 - 500x Digital Portable Metallurgical Optical Microscope A31.2501



Small & Simple, with bright white LED light built-in, and high quality image!

A13.2501 Portable Metallurgical Microscope		
Total Magnification	100-400x	
Head	Monocular	
Eyeiece	WF10X	
LWD Plan Objective (no cover glass)	PLL 10X/0.25 WD:7.3mm	
	PLL 40X/0.65 WD:0.5mm	
Mechanical Tube Length	160mm	
Focusing	Coaxial coarse/fine focus system, range:36mm	
Light Source	Coaxial illumination, build-in LED rechargeable, brightness adjustable	
Packing	Portable aluminium case	
A13.2501 Portable Metallurgical Microscope Optional Accessories		
Adapter	For Canon SLR Digital Camera	A55.2502
Stand	Magnet Stand	A54.2501
LWD Plan Objective (no cover glass)	PLL 20x	A5M.0234-20
	PLL 60x	A5M.0234-60
	PLL 80x	A5M.0234-80
	PLL 100x	A5M.0234-100



A13.2501





Opto-Edu (Beijing) Co., Ltd.

0086 13911110627    sale@optoedu.com    cnoec.com  
F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China