



Infinity Trinocular Biological Microscope 6V 30W Halogen Lamp A12.0204

Basic Information

- Place of Origin: China
- Brand Name: CNOEC
- Certification: CE
- Model Number: A12.0204
- Minimum Order Quantity: 1
- Price: Negotiation
- Packaging Details: Carton
- Delivery Time: 1 month
- Payment Terms: T/T, West Union and Paypal
- Supply Ability: 10000




Product Specification

- Head: Binocular Or Trinocular Head Inclined 30° Available
- Light: Kohler Illumination With 6V 20W Halogen Lamp
- Highlight: high powered microscope, binocular biological microscope



More Images



Product Description

Infinity Trinocular Biological Microscope 6V 30W Halogen Lamp A12.0204

New design in 2009!

- Infinity optical system!
- Binocular or trinocular head inclined 30° available
- Coaxial coarse/fine focusing system with tension adjustable and up stop
- Kohler illumination with 6V 20W halogen lamp



OPTO-EDU (BEIJING) CO., LTD.

F-1501 Wanda Plaza, No.18 Shijingshan Road, Beijing 100043, China
Tel:+8610 88696020 Fax:+8610 88696085

A12.0204 Laboratory Biological Microscope



Product Details



A5C.0210-1 Small Turret
Phase Contrast Condenser



A5C.0210-4 Big Turret
Phase Contrast Condenser



A5C.0210-2 Flash Board
Phase Condenser Set



A5C.0210-3 Pull Board
Phase Condenser Set



A12.0204 Trinocular Laboratory Microscope



A5C.0234 Infinity Plan Phase Contrast Objective
A5C.0212-1 Centering Telescope 10x



A5D.0210-D Dry Dark Field Condenser
A5D.0210-4 Immersion Dark Field Condenser




Binocular Head, WF10x/22mm High
Eyepoint Eyepiece, Diopter Adjustable




Infinity Plan Achromatic Objective 4x/0.1,
10x/0.25, 40x/0.65(S), 100x/1.25(O.S)



Light Collector For Halogen Lamp, With
Field Diaphragm



A12.0204 Laboratory Microscope System Size



A12.0204 Laboratory Microscope System Diaphragm

	A12.0204-A	A12.0204-B
Optical System	Infinity Optical System	
Head	Trinocular Head, 30° Inclined, Interpupillary Distance 55~75mm, Both Eyepiece Tube Diopter Adjustable +/-5°, 100% Light Pass To Photo Tube 360° Rotatable Head Adaptor For Easy Viewing Height Adjustable	
Eyepiece	WF10x/22mm Wide Field Eyepiece	
Nosepiece	Backward Quadruple Nosepiece, Ball Bearing Inner Locating	Backward Quintuple Nosepiece, Ball Bearing Inner Locating
Objective	Infinity Plan Achromatic Objective: PL4x/0.10, W.D.21mm PL10x/0.25, W.D. 5.0mm PL40x/0.65, W.D. 0.66mm (S) PL100x/1.25, W.D.0.36mm (S,O)	Infinity Plan Achromatic Objective: PL4x/0.10, W.D.21mm PL10x/0.25, W.D. 5.0mm PL20x/0.40, W.D.8.8mm PL40x/0.65, W.D. 0.66mm (S) PL100x/1.25, W.D.0.36mm (S,O)
Working Stage	Double Layer Mechanical Stage 210*140mm, Moving Range 75*50mm	
Focusing	Coaxial Coarse & Fine Focus Knob, With Tension Adjustable And Limit Stopper, Minimum Division 0.002mm	
Condenser	Abbe Condenser N.A.1.25, Rack & Pinion Height Adjustable	
Filter	Blue Filter, Ground Glass	
Light Source	Kohler Illumination, 6V30W Halogen Lamp, Brightness Adjustable, Light Collector With Iris Diaphragm	

A12.0204 Laboratory Biological Microscope Optional Accessories		
Eyepiece	Retical Eyepiece, WF10x/22mm, Micrometer Div 0.1mm	A51.0205-G10B
Objective	Infinity Plan Achromatic Objective PL60x/0.85, W.D.0.25mm	A52.0206-60
Filter	Green Filter	A56.0209-02
	Yellow Filter	A56.0209-03
Phase Contrast	Centering Telescope	A5C.0212-2
	Infinity Plan Phase Contrast Objective 10x/0.25 PHP	A5C.0234-P10
	Infinity Plan Phase Contrast Objective 20x/0.40 PHP, Spring	A5C.0234-P20
	Infinity Plan Phase Contrast Objective 40x/0.65 PHP, Spring	A5C.0234-P40
	Infinity Plan Phase Contrast Objective 100x/1.25 PHP, Spring, Oil	A5C.0234-P100
	Small Turret Phase Contrast Condenser	A5C.0210-1
	Big Turret Phase Contrast Condenser	A5C.0210-4
	Flashboard Phase Contrast Condenser With Slides *4	A5C.0210-2
	Pullboard Phase Contrast Condenser With Sildes * 2	A5C.0210-3
Polarizing	Simple Polarizing Set Polarizer, 360° Rotatable, Used On Base Light Collector Analyzer Slide, Inserted Into Hole On Nosepiece	A5P.0201-2
Dark Field	Dark Field Condenser, Dry	A5D.0210-D
	Dark Field Condenser, Immersion	A5D.0210-I
Adapter	CCD Adapter 0.4x	A55.0202-01
	CCD Adapter 0.5x	A55.0202-04
	CCD Adapter 1.0x	A55.0202-02
	CCD Adapter 0.5x With Dividing 0.1mm/Divv	A55.0202-03
	Digital SLR Camera Adapter for Canon EF Series	A55.0204-02





A12.0204 Laboratory Microscope



A5C.0210-1 Small Turret
Phase Contrast Condenser



A5C.0210-4 Big Turret
Phase Contrast Condenser



A5C.0210-2 Flash Board
Phase Condenser Set



A5C.0210-3 Pull Board
Phase Condenser Set



A5C.0234 Infinity
Plan Phase Contrast Objective
A5C.0212-1
Centering Telescope 10x



A5D.0210-D
Dry Dark Field Condenser
A5D.0210-I
Immersion Dark Field Condenser





Binocular Head, WF10x/22mm
High Eyepoint Eyepiece,
Diopter Adjustable




Infinity Plan Achromatic
Objective 4x/0.1, 10x/0.25,
40x/0.65(S), 100x/1.25(O,S)




Light Collector For Halogen
Lamp, With Field Diaphragm



A12.0204 Laboratory Microscope System Size





Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China