



Binocular 40X - 1000X Halogen 6V/20W Biological Microscope A12.1106

Basic Information

- Place of Origin: China
- Brand Name: CNOEC
- Certification: CE
- Model Number: A12.1106
- Minimum Order Quantity: 1
- Price: Negotiation
- Packaging Details: Carton
- Delivery Time: 1 month
- Payment Terms: T/T, West Union and Paypal
- Supply Ability: 10000



Product Specification

- Highlight: binocular biological microscope, biological compound microscope



Product Description

Binocular 40X - 1000X Halogen 6V/20W Biological Microscope A12.1106

A12.1106 Biological Microscope		-A	-B
Head	Compensation Free Binocular,30°inclined,55mm-75mm	●	-
	Compensation Free Trinocular,30°inclined,55mm-75mm	-	●
Eyeiece	WF10X	●	●
Objective	Achromatic 4X	●	●
	Achromatic 10X	●	●
	Achromatic 40X, Spring	●	●
	Achromatic 100X, Spring,Oil	●	●
Nosepiece	Quadruple Revolving Nosepiece	●	●
Stage	Double Layer Mechanical Stage,Size140x140mm,Moving Range:75mmx50mm	●	●
Focusing	Coaxial Coarse and Fine Focusing with Rack and Pinion Mechanism; Fine Focusing Scale Value 0.002mm	●	●
Condenser	ABBE N.A.1.25Condenser with Iris Diaphragm & Filter	●	●
Illumination	High Brightness Kohler Illumination	●	●
Light Source	Halogen Bulb 6V/20W AC85V-230V,,Adjustble Brightness	●	●
Optional Accessories			
Eyeiece	WF16X		
	WF20X		
	P16X		
Objective	Achromatic 20X		
	Achromatic 60X, Spring		
	Semi-plan Objective 4X		
	Semi-plan Objective 10X		
	Semi-plan Objective 20X		
	Semi-plan Objective 40X(S)		
	Semi-plan Objective 60X(S)		
	Semi-plan Objective 100X(S)		
	Plan Objective 4X		
	Plan Objective 10X		
	Plan Objective 20X		
	Plan Objective 40X(S)		
	Plan Objective 60X(S)		
	Plan Objective 100X(S)		
	Infinity Plan Objective 4X		
	Infinity Plan Objective 10X		
	Infinity Plan Objective 20X		
	Infinity Plan Objective 40X(S)		
	Infinity Plan Objective 60X(S)		
	Infinity Plan Objective 100X(S)		
Attachment	Phase Contrast Attachment		
Dark Field	Darkfield		
Photo Adapter	Photo Adapter		
CCD Adapter	0.5x		
	0.57x		
	0.75x		



A12.1106-A



A12.1106-B



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China