



Infinity Trinocular Biological Microscope Kohler Illumination A12.1107

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC
- Certification: CE
- Model Number: A12.1107
- Minimum Order Quantity: 1
- Price: Negotiation
- Packaging Details: Carton
- Delivery Time: 1 month
- Payment Terms: T/T, West Union and Paypal
- Supply Ability: 10000



Product Specification

- Highlight: high powered microscope, biological compound microscope



More Images



for more products please visit us on cnoec.com

Product Description

Infinity Trinocular Biological Microscope Kohler Illumination A12.1107

A12.1107 Biological Microscope		-A	-B	-C	-D
Head	Compensation Free Binocular,30°inclined,55mm-75mm	●	-	●	-
	Compensation Free Trinocular,30°inclined,55mm-75mm	-	●	-	●
Eyepiece	WF10X /18mm	●	●	-	-
	WF10X /22mm	-	-	●	●
Objective	Achromatic 4X	●	●	-	-
	Achromatic 10X	●	●	-	-
	Achromatic 40X, Spring	●	●	-	-
	Achromatic 100X, Spring,Oil	●	●	-	-
	Infinity Objective 4X	-	-	●	●
	Infinity Objective 10X	-	-	●	●
	Infinity Objective 40X, Spring	-	-	●	●
	Infinity Objective 100X, Spring,Oil	-	-	●	●
Nosepiece	Quadruple Revolving Nosepiece	●	●	●	●
Stage	Double Layer Mechanical	●	●	●	●
	Stage,Size180mmx150mm,Moving Range:75mmx50mm	●	●	●	●
Focusing	Coaxial Coarse and Fine Focusing with Rack and Pinion Mechanism; Fine Focusing Scale Value 0.002mm	●	●	●	●
Condenser	ABBE N.A.1.25Condenser with Iris Diaphragm & Filter	●	●	●	●
Illumination	High Brightness Kohler Illumination	●	●	●	●
Light Source	Halogen Bulb 6V/20W,AC85V-230V, Adjustable Brightness	●	●	●	●
Optional Accessories					
Eyepiece	WF16X				
	WF20X				
	P16X				
Objective	Achromatic 20X				
	Achromatic 60X, Spring				
	Semi-plan Objective 4X				
	Semi-plan Objective 10X				
	Semi-plan Objective 20X				
	Semi-plan Objective 40X(S)				
	Semi-plan Objective 60X(S)				
	Semi-plan Objective 100X(S)				
	Plan Objective 4X				
	Plan Objective 10X				
	Plan Objective 20X				
	Plan Objective 40X(S)				
	Plan Objective 60X(S)				
	Plan Objective 100X(S)				
Condenser	Swing ABBE N.A.1.25Condenser with Iris Diaphragm & Filter				
Dark Field	Darkfield				
Photo Adapter	Photo Adapter				
CCD Adapter	0.5x				
	0.57x				
	0.75x				



A12.1107-A



A12.1107-B



Back side of microscope body



Compensation Free
Binocular Head



ABBE N.A. 1.25 Condenser
with Iris Diaphragm & Filter





A12.1107-C



A12.1107-D



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China