



Laboratory Research Trinocular Polarizing Light Microscope With CE A15.1103

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A15.1103
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Highlight: **student compound microscope ,
highest magnification microscope**

Product Description

Laboratory Research Trinocular Polarizing Light Microscope With CE A15.1103

| | | | | | |
|--|---|-----|---|---|----|
| A15.1103 Polarizing Microscope Specification | | | - | A | -B |
| Head | Compensation Free Binocular,30°inclined,55mm-75mm | | ● | - | |
| | Compensation Free Trinocular,30°inclined,55mm-75mm | | - | ● | |
| Eyepiece | Wide field WF10X/22mm | | ● | ● | |
| | Gridding Eyepiece | | ● | ● | |
| | Cross Eyepiece | | ● | ● | |
| | Eyepiece with Reticule | | ● | ● | |
| Objective | Stress-Free Infinity Plan-Achromatic Metallurgical Objective | 40X | | ● | ● |
| | | 60X | | ● | ● |
| | Stress-Free Infinity Plan Achromatic Objective | 4X | ● | ● | |
| | | 10X | ● | ● | |
| | | 20X | ● | ● | |
| | | 40X | ● | ● | |
| | | 60X | ● | ● | |
| Condenser | Swing Type N.A.1.25 Abbe Condenser with Iris Diaphragm | | ● | ● | |
| Focus System | Coaxial Coarse and Fine Focusing with Rack and Pinion Mechanism; Fine Focusing Scale Value 0.002mm. | | ● | ● | |
| Nosepiece | Quintuple Nosepiece with centering adjustable | | ● | ● | |
| Filter | Green Inserting Type Filting | | ● | ● | |
| | Blue Inserting Type Filting | | ● | ● | |
| | Neutral Inserting Type Filting | | ● | ● | |
| Stage | Rotary Stage,Stage Size: Φ172mm,Moving Range:360° | | ● | ● | |
| Testing Piece | Gypsum 1 λTesting Piece | | ● | ● | |
| | Mica1/4 λ | | ● | ● | |
| | Quartz Wedge Testing Piece | | ● | ● | |
| Illumination | High brightness Fixed Illumination | | ● | ● | |
| Light Source | Upper:Halogen Bulb 12V/30W,AC 85V-230V ,Adjustble Brightness | | ● | ● | |
| | Down:Halogen Bulb 12V/30W,AC 85V-230V ,Adjustble Brightness | | ● | ● | |
| Checking Tool | 0.01mm Micrometer | | ● | ● | |
| Optional Accessories | | | | | |
| Eyepiece | WF16X | | | | |
| | WF20X | | | | |
| | P16X | | | | |
| Stress-Free Infinity Plan-Achromatic Metallurgical Objective | 4X | | | | |
| | 10X | | | | |
| | 20X | | | | |
| Photo Adapter | Photo Adapter | | | | |
| CCD Adapter | 0.5X | | | | |
| | 0.57X | | | | |
| | 0.75X | | | | |



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China