



Inclined Head A23.1001 Stereo Zoom Binocular Microscope Max 0.8~8x

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: OPTO-EDU
- Certification: CE, Rohs
- Model Number: A23.1001
- Minimum Order Quantity: 1pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: Western Union, T/T Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Head: Binocular Head 20° Inclined
- Eyepiece: Extra Wide Field Eyepiece
EWF10x/22mm Diopter
Adjustable, Interpupillary 55-75mm
- Auxiliary Lens: Plan Achromatic Lens 1.0x
- Focusing: Coarse Focusing Konb, Range 105mm
- Illumination: Transmit & Reflect LED Illumination,
Brightness Adjustable, Working Life 6000
Hours
- Working Distance: 78mm
- Highlight: **0.8x stereo zoom binocular microscope,
inclined head stereo zoom binocular
microscope
, 8x zoom stereo microscope**

Product Description

Parallel Optical System, Module Design to Add Function Easily
Large Zoom Ratio Up to 1:10 With Tilting Head 0°-35° Viewing Angle Adjustable for User
Transmit & Reflect LED Illumination, Brightness Adjustable, Working Life 6000 Hours
Complete Accessories Include Polarizing, Dark Field, Fluorescent Attachment
Plan Apochromatic Objective 0.5x Optional



OPTO-EDU (BEIJING) CO., LTD.

OPTO-EDU

F-1501 Wanda Plaza, No.18 Shijingshan Road, Beijing 100043, China
Tel:+8610 88696020 Fax:+8610 88696085

A23.1001 Series

Zoom Stereo Microscope, Max 0.8~8x



Product Details



Product Details





A23.1001 Standard Set

**A23.1001 Standard Set with
A55.1003-S1 Photo Tube, Beam Splitter 1 Port**

Product Details



A5F.1020 EPI-Fluorescent Attachment For Stereo Microscope



A23.1001
Zoom Stereo Microscope, Max 0.8~8x



A23.1001-F
Flourescent Zoom Stereo Microscope


Product Details



A23.1001-T Standard Set

**A23.1001-T Standard Set with
A55.1003-S1 Photo Tube, Beam Splitter 1 Port**

Product Dimension




A23.1001

A23.1001-T

A23.1001-F

Accessories



Product Specification

A23.1001 Zoom Stereo Microscope		A23.1 001- 06	A23.1 001- 08	A23.1 001- 10
Optical System	Parallel Optical System, Module Design to Add Function Easily	●	●	●
Head	Binocular Head 20° Inclined	●	●	●
	Titling Head 0°-35° Viewing Angle Adjustable for User With Titling Head, Model Change to A23.1001-T	o	o	o
Eyepiece	Extra Wide Field Eyepiece EWF10x/22mm,Diopter Adjustable,Interpupillary 55-75mm	●	●	●
Zoom Body	0.8~5.0x, Zoom Ratio 1:6	●	-	-
	0.8~6.4x, Zoom Ratio 1:8	-	●	-
	0.8~8.0x, Zoom Ratio 1:10	-	-	●
Auxiliary Lens	Plan Achromatic Lens 1.0x	●	●	●
Working Distance	78mm	●	●	●
Focusing	Coarse Focusing Konb, Range 105mm	●	●	
	Coaxial Coarse & Fine Focusing Knob, Fine Division 0.015mm, Coarse Stroke 38.4mm/Rotation, Fine Stroke 1.5mm/Rotation, Moving Range 105mm , Focusing Mount Inner Dia. 76mm			●
Illumination	Transmit & Reflect LED Illumination, Brightness Adjustable, Working Life 6000 Hours			
A23.1001 Zoom Stereo Microscope		Cata. No.		
Head	Titling Head 0°-35° Viewing Angle Adjustable for User Model Change to A23.1005	A53.1014		
Eyepiece	WF15x/16mm	A51.1027-1516		
	WF20x/12mm	A51.1027-2012		
	WF30x	A51.1027-30		
Auxiliary Lens	Plan Apochromatic Lens 0.5x	A52.1023-05A		
	Plan Achromatic Lens 1.0x	A52.1023-10		
	Plan Achromatic Lens 2.0x	A52.1023-20		
	Achromatic Lens 0.3x	A52.1023-03		
	Achromatic Lens 0.5x	A52.1023-05		
Illumination	Coaxial Illumination Set:	A56.1021		
	--Cold Light Power Supplier,			
	--Dual Light Guide,			
	--Coaxial Episcopic Illuminator,			
	--Iris Diaphragm	A56.1021-2		
	Iris Diaphragm			
Working Stage	Temperature Control Stage	A54.1004		
	Movable Working Stage	A54.1007		
Stand Click to See	Universal Table Stand	A54.1011		
		A54.1012		
	Track Stand, Transmit & Reflect Light Source With Coarse & Fine Focusing Optional, Model -F, Like A54.1019-F	A54.1013		
		A54.1014		
		A54.1019		
		A54.1020		
Fluorescent	EPI-Fluorescent Attachment For Stereo Microscope Model Change to A36.100-F	A5F.1020		
Dark Field	Dark Field Attachment For Stereo Microscope	A5D.1001		
Polarizing	Simple Polarizing Attachment For Stereo Microscope	A5P.1003		
LCD Monitor	8.4" TFT Touch Screen LCD Monitor & Connector Model Change to A36.1001	A59.1006		
Adapter	Photo Tube, Beamsplitter 1 Port	A55.1003-S1		
	Photo Tube, Beamsplitter 2 Port	A55.1003-S2		
	Photo Adapter, 2.0x, For Nikon or Canon Dlgital SLR Camera	A55.1020-2X		
	C-mount 0.5x	A55.1021-05		
	C-mount 1.0x	A55.1021-10		



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China