



A23.3667L OPTO-EDU 6.7x - 45x Stereoscopic Zoom Microscope With Pole Stand

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: OPTO-EDU
- Certification: CE, Rohs
- Model Number: A23.3667L
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Eyepiece: WF10x/20mm
- Zoom Lens: 0.67~4.5x
- Zoom Ratio: 1:6.7
- Magnification: 6.7x-45x
- Working Distance: 100mm, With Auxiliary Objectives 26~287mm
- Stand: Pole Stand, No Light Source
- Highlight: 45x stereoscopic zoom microscope, pole stand stereoscopic zoom microscope, 6.7x opto edu microscope



More Images



Product Description

Eyepiece WF10X/20mm, High Eyepoint, Eyepiece Tube Diopter Adjustable
0.67x-4.5x, Magnification Up To 3.35x-270x With Optional Eyepieces And Auxiliary Objectives
A54.3640 35 Kinds Stands Optional For Different Usage & Condition

A23.3667L-B1 Zoom Stereo Microscope Specification		B	T
Head	Binocular Head, 45° Inclined, 360°Rotatable, ±5 Diopter Adjustable, Interpupillary Adjustment 54-75mm	O	
	Trinocular Head, 45° Inclined, 360°Rotatable, ±5 Diopter Adjustable, Interpupillary Adjustment 54-75mm, Eyepiece Tube Adapter		O
Eyepiece	WF10x/20mm,High Eyepoint		
Zoom Range	0.67~4.5x		
Zoom Ratio	1:6.7		
Magnification	6.7~45x, (2x~225x With Auxliary Objective)		
Working Distance	100mm, With Auxiliary Objectives 26~287mm		
Stand	Pole Stand, No Light Source		
Optional Accessories		Cata.No.	
Head	With WF10X eyepiece; 6.7-45x continuously adjustable	A53.3667L-B	
	Trinocular With WF10X eyepiece; 6.7-45x continuously adjustable	A53.3667L-T	
Eyepiece	WF10x/20mm, High-eyepoint	A51.3621-1020	
	WF10x/22mm, High-eyepoint	A51.3621-1022	
	WF10x/23mm, High-eyepoint	A51.3621-1023	
	WF15x/15mm, High-eyepoint	A51.3621-1515	
	WF20x/10mm, High-eyepoint	A51.3621-2010	
	WF25x/9mm, High-eyepoint	A51.3621-2509	
	WF30x/8mm, High-eyepoint	A51.3621-3008	
	WF10x/22mm, High-eyepoint, Diopter Adjustable	A51.3622-1022	
	WF15x/15mm, High-eyepoint, Diopter Adjustable	A51.3622-1515	
	WF20x/10mm, High-eyepoint, Diopter Adjustable	A51.3622-2010	
	WF25x/9mm, High-eyepoint, Diopter Adjustable	A51.3622-2509	
	WF30x/8mm, High-eyepoint, Diopter Adjustable	A51.3622-3008	
Auxiliary Lens	0.3X,W.D. 287mm	A52.3622-0.3x	
	0.5X,W.D. 177mm	A52.3622-0.5x	
	0.7X,W.D. 113mm	A52.3622-0.7x	
	0.75X,W.D. 116mm	A52.3622-0.75x	
	1.5X,W.D. 47mm	A52.3622-1.5x	
	2.0X,W.D. 26mm	A52.3622-2.0x	
CCD Adapter	1.0x, For CTV	A55.3622-1.0x	
	0.50x, For 1/2" CCD	A55.3622-0.5x	

A23.3667L Details



Eyepiece WF10X/20mm, High Eyepoint, Eyepiece Tube Diopter Adjustable



0.67x-4.5x, Magnification Up To 3.35x-270x With Optional Eyepieces And Auxiliary Objectives



A23.3667L Series Binocular Models



Binocular Head



35 Stands



A23.3667L-B1



A23.3667L-B2



A23.3667L-B2L



A23.3667L-B3



A23.3667L-B3LS



A23.3667L-B2LD



A23.3667L-B4



A23.3667L-B5



A23.3667L-B6



A23.3667L-B7



A23.3667L-B8L



A23.3667L-B8LD



A23.3667L-B9L



A23.3667L-B9LD



A23.3667L-B8LS



A23.3667L-L



A23.3667L-2L



A23.3667L-B9LS



A23.3667L-J1



A23.3667L-J2



A23.3667L-J3L



A23.3667L-J4L



A23.3667L-J4LD



A23.3667L-J3LD



A23.3667L-STL1



A23.3667L-STL2



A23.3667L-STL3



A23.3667L-STL4



A23.3667L-STL5



A23.3667L-STL8



A23.3667L-STL6A



A23.3667L-STL6B



A23.3667L-STL6C

A23.3667L Series Trinocular Models



Trinocular Head



35 Stands



A23.3667L-B1T



A23.3667L-B2T



A23.3667L-B2LT



A23.3667L-B3T



A23.3667L-B3LST



A23.3667L-B2LDT



A23.3667L-B4T



A23.3667L-B5T



A23.3667L-B6T



A23.3667L-B7T



A23.3667L-B8LT



A23.3667L-B8LDT



A23.3667L-B9LT



A23.3667L-B9LDT



A23.3667L-B8LST



A23.3667L-LT



A23.3667L-2LT



A23.3667L-B9LST



A23.3667L-J1T



A23.3667L-J2T



A23.3667L-J3LT



A23.3667L-J4LT



A23.3667L-J4LDT



A23.3667L-J3LDT



A23.3667L Zoom Stereo Microscope Optical Data Sheet							
Auxiliary Lens	0.3X	0.4X	0.5X	0.75X	1.0X	1.5X	2.0X
W.D.	287mm	217mm	177mm	117mm	100mm	47mm	26mm
WF10X	Magnification	2.0-13.5x	2.7-18.0x	3.4-22.5x	5.0-33.8x	6.7-45.0x	10.1-67.5x
	View Field(mm)	109.5-16.3	82.1-12.2	65.7-9.8	43.8-6.5	32.8-4.9	21.9-3.3
WF15X	Magnification	3.0-20.3x	4.0-27.0x	5.0-33.8x	7.5-50.6x	10.1-67.5x	15.1-101.3x
	View Field(mm)	79.6-11.9	59.7-8.9	47.8-7.1	31.8-4.7	23.9-3.6	15.9-2.4
WF20X	Magnification	4.0-27.0x	5.4-36.0x	6.7-45.0x	10.1-67.5x	13.4-90.0x	20.1-135.0x
	View Field(mm)	59.7-8.9	44.8-6.7	35.8-5.3	12.9-2.4	17.9-2.7	11.9-1.8
WF25X	Magnification	5.0-33.8x	6.7-45.0x	8.4-56.3x	12.6-84.4x	16.8-112.5x	25.1-168.8x
	View Field(mm)	44.8-6.7	33.6-5.0	26.9-4.0	17.9-2.7	13.4-2.0	9.0-1.3



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China