



A23.3645N OPTO-EDU Stereo Optical Microscope 0.7~4.5x Trinocular 1/6.5

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: OPTO-EDU
- Certification: CE, Rohs
- Model Number: A23.3645N
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Eyepiece: WF10x/20mm
- Zoom Lens: 0.7~4.5x
- Zoom Ratio: 1:6.5
- Magnification: 7x-45x, Up To 3.5x-180x With Optional Eyepieces And Auxiliary Objectives
- Working Distance: 100mm, With Auxiliary Objectives 30~165mm
- Field Of View: 1~95.2mm
- Stand: Pole Stand, No Light Source
- Highlight: trinocular stereo optical microscope, 4.5x stereo optical microscope, opto edu stereo zoom microscopes



More Images



Product Description


Eyepiece WF10X/20mm, High Eyepoint, Eyepiece Tube Diopter Adjustable
 0.7x-4.5x, Magnification Up To 3.5x-180x With Optional Eyepieces And Auxiliary Objectives
 A54.3640 35 Kinds Stands Optional For Different Usage & Condition
 Hinged triad, 50/50%, 20/80%, WF10X/20mm High Eyepieces, With A Variety Of Base Options, Up to 3.5-270x With Optional Eyepieces And Auxiliary Objectives.

A23.3645N-B1 Zoom Stereo Microscope Specification		B	T
Head	Binocular Head, 45° Inclined, 360°Rotatable, ±5 Diopter Adjustable, Interpupillary Adjustment 54-75mm	O	
	Trinocular Head, 45° Inclined, 360°Rotatable, ±5 Diopter Adjustable, Interpupillary Adjustment 54-75mm, Eyepiece Tube Adapter		O
Eyepiece	WF10x/20mm, Both With Eye-cup		
Zoom Lens	0.7~4.5x		
Zoom Ratio	1:6.5		
Magnification	7x-45x, Up To 3.5x-180x With Optional Eyepieces And Auxiliary Lens		
Working Distance	100mm, With Standard 1x Lens		
Field Of View	1~95.2mm		
Stand	Pole Stand, No Light Source		
Optional Accessories		Cata.No.	
Head	Binocular Head + WF10x/20mm	A53.3645N-B	
	Trinocular Head + WF10x/20mm	A53.3645N-T	
Eyepiece	WF10x/20mm, High-eyepoint	A51.3621-1020	
	WF10x/22mm, High-eyepoint	A51.3621-1022	
	WF10x/23mm, High-eyepoint	A51.3621-1023	
	WF15x/15mm, High-eyepoint	A51.3621-1515	
	WF20x/10mm, High-eyepoint	A51.3621-2010	
	WF25x/9mm, High-eyepoint	A51.3621-2509	
	WF30x/8mm, High-eyepoint	A51.3621-3008	
	WF10x/22mm, High-eyepoint, Diopter Adjustable	A51.3622-1022	
	WF15x/16mm, High-eyepoint, Diopter Adjustable	A51.3622-1516	
	WF20x/12mm, High-eyepoint, Diopter Adjustable	A51.3622-2012	
	WF25x/9mm, High-eyepoint, Diopter Adjustable	A51.3622-2509	
	WF30x/8mm, High-eyepoint, Diopter Adjustable	A51.3622-3008	
Auxiliary Lens	0.5X,W.D. 165mm	A52.3621-0.5x	
	1.5X,W.D. 45mm	A52.3621-1.5x	
	2.0X,W.D. 30mm	A52.3621-2.0x	
CCD Adapter	1.0x, For CTV	A55.3621-1.0x	
	0.50x, For 1/2" CCD	A55.3621-0.5x	
	0.35x, For 1/3" CCD	A55.3621-0.35x	
	0.50x, For 1/2" CCD	A55.3621-N0.5x	

A23.3645N Details



Eyepiece WF10X/20mm, High Eyepoint, Eyepiece Tube Diopter Adjustable



0.7x-4.5x, Magnification Up To 3.5x-180x With Optional Eyepieces And Auxiliary Objectives



A23.3645N Series Binocular Models



Binocular Head



35 Stands



A23.3645N-B1



A23.3645N-B2



A23.3645N-B2L



A23.3645N-B3



A23.3645N-B3LS



A23.3645N-B2LD



A23.3645N-B4



A23.3645N-B5



A23.3645N-B6



A23.3645N-B7



A23.3645N-B8L



A23.3645N-B8LD



A23.3645N-B9L



A23.3645N-B9LD



A23.3645N-B8LS



A23.3645N-L



A23.3645N-2L



A23.3645N-B9LS



A23.3645N-J1



A23.3645N-J2



A23.3645N-J3L



A23.3645N-J4L



A23.3645N-J4LD



A23.3645N-J3LD



A23.3645N-STL1



A23.3645N-STL2



A23.3645N-STL3



A23.3645N-STL4



A23.3645N-STL5



A23.3645N-STL8



A23.3645N-STL6A



A23.3645N-STL6B



A23.3645N-STL6C

A23.3645N Series Trinocular Models



Trinocular Head



35 Stands



A23.3645N-B1T



A23.3645N-B2T



A23.3645N-B2LT



A23.3645N-B3T



A23.3645N-B3LST



A23.3645N-B2LDT



A23.3645N-B4T



A23.3645N-B5T



A23.3645N-B6T



A23.3645N-B7T



A23.3645N-B8LT



A23.3645N-B8LDT



A23.3645N-B9LT



A23.3645N-B9LDT



A23.3645N-B8LST



A23.3645N-LT



A23.3645N-2LT



A23.3645N-B9LST



A23.3645N-J1T



A23.3645N-J2T



A23.3645N-J3LT



A23.3645N-J4LT



A23.3645N-J4LDT



A23.3645N-J3LDT



A23.3645N-STL1T



A23.3645N-STL2T



A23.3645N-STL3T



A23.3645N-STL4T



A23.3645N-STL5T



A23.3645N-STL8T



A23.3645N-STL6AT



A23.3645N-STL6BT






A23.3645N-STL6CT

Technical Parameters:

A23.3645N Zoom Stereo Microscope Optical Data Sheet									
Objective	Standard 1.0x		Auxiliary 0.5X		Auxiliary 1.5X		Auxiliary 2X		
W.D.	100mm		165mm		45mm		30mm		
Eyepiece	Mag	View Field	Mag	View Field	Mag	View Field	Mag	View Field	
	7.0X	28.6mm	3.5X	57.1mm	10.5X	19mm	14X	14.3mm	
10X/20mm	45.0X	4.4mm	22.5X	8.9mm	67.5X	3mm	90X	2.2mm	
	10.5X	21.1mm	5.25X	42.8mm	15.75X	14.3mm	21X	10.7mm	
15X/15mm	67.5X	3.3mm	33.75X	6.7mm	101.25X	2.2mm	135X	1.7mm	
	14.0X	14.3mm	7.0X	28.6mm	21.0X	9.5mm	28X	7.1mm	
20X/10mm	90.0X	2.2mm	45.0X	4.4mm	135.0X	1.5mm	180X	1.1mm	



 0086 13911110627  sale@optoedu.com  cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China