



A22.3660 OPTO EDU Stereo Microscope Binocular 1x3x Or 2X/4X

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: OPTO-EDU
- Certification: CE, Rohs
- Model Number: A22.3660
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month

OPTO-EDU



Product Specification

- Eyepiece: WF10x/20mm
- Head: 45°Binocular Head
- Objective: Turret Objective 1x3x Or 2X/4X
- Magnification Range: 10X-80X
- Working Distance: 100mm
- Stand: Pole Stand, No Light Source
- Highlight: wf10x stereo microscope binocular, 20mm stereo microscope binocular, stereo binocular opto edu microscope



More Images



Product Description

Stereo microscope is widely used in materials such as macroscopic surface observation, failure analysis, fracture analysis and other industries field

Binocular Head, 45° Inclined, 360° Rotatable

Wide Field, High-Eyepoint Eyepiece-WF10X20mm

Turret Objective 1x/3x, 2x/4x or 1x/2x

Pole Stand, No Light Source

A22.3660-B1 Stereo Microscope Specification	
Head	45° Binocular head
Eyepiece	Wide Field, High-eyepoint Eyepiece-WF10X20mm
Objective	Turret Objective 1x3x or 2X/4X
Magnification Range	10X-80X
Inter-Pupillary Distance	55-75mm
Diopter Range	5mm
Working Distance	100mm
Stand	Pole Stand, No Light Source

A22.3660 Details



Eyepiece WF10X/20mm, High Eyepoint, Eyepiece Tube Diopter Adjustable



20x/40x, Magnification Up To 10x.80x With Optional Eyepieces And 1X2X, 1X3X, 2X4X, Objectives

A22.3660 Series Models



Binocular Head



35 Stands



A22.3660-B1



A22.3660-B2



A22.3660-B2L



A22.3660-B3



A22.3660-B3LS



A22.3660-B2LD



A22.3660-B4



A22.3660-B5



A22.3660-B6



A22.3660-B7



A22.3660-B8L



A22.3660-B8LD



A22.3660-B9L



A22.3660-B9LD



A22.3660-B8LS



A22.3660-L



A22.3660-2L



A22.3660-B9LS



A22.3660-J1



A22.3660-J2



A22.3660-J3L



A22.3660-J4L



A22.3660-J4LD



A22.3660-J3LD



A22.3660-STL1



A22.3660-STL2



A22.3660-STL3



A22.3660-STL4



A22.3660-STL5



A22.3660-STL8



A22.3660-STL6A




A22.3660-STL6B



A22.3660-STL6C

A22.36, A23.36 Total 8 Stereo Microscope Series



Objective Optics Parameter

Model	Magnification	Eyepiece			Working Distance
		SZMEWh 10X20	SZMEWh 15X15	SZMEWh 20X10	
A22.3660-B1-12	1X	10X20mm	15X15mm	20X10mm	100mm
	2X	10X20mm	30X7.5mm	40X5mm	
A22.3660-B1-13	1X	10X20mm	15X15mm	20X10mm	
	3X	30X6.7mm	45X5mm	60X3.3mm	
A22.3660-B1-24	2X	10X20mm	30X7.5mm	40X5mm	
	4X	40X5mm	60X3.75mm	80X2.5mm	



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China