



Stereo Optical Microscope Parellel Zoom Stereo 0.8 - 8x , 0.8 - 6.4x A23.1001

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC
- Certification: CE, Rohs
- Model Number: A23.1001
- Minimum Order Quantity: 1
- Price: \$500-1000
- Packaging Details: Carton Packing
- Delivery Time: 5-30
- Payment Terms: L/C, T/T, Western Union
- Supply Ability: 10000



Product Specification

- Highlight: stereo zoom microscopes, stereo binocular microscopes



More Images



Product Description

Stereo Optical Microscope Parellel Zoom Stereo 0.8 - 8x , 0.8 - 6.4x A23.1001

A23.1001 Parallel Zoom Stereo Microscope, the newest advanced zoom stereo microscope, can be used extensively for scientific research, industry, biological and medical field.

- Parallel Optical System, Module Design to Add Function Easily
- Large Zoom Ratio Up to 1:10
- Transmit & Reflect LED Illumination, Brightness Adjustable, Working Life 6000 Hours
- Complete Accessories Include Polarizing, Dark Field, Fluorescent Attachment
- Plan **A**pochromatic Objective 0.5x Optional

Our Product Introduction

A23.1001 Zoom Stereo Microscope		A	B	C
Optical System	Parallel Optical System, Module Design to Add Function Easily			
Head	Binocular Head 20° Inclined			
Eyeiece	Extra Wide Field Eyeiece EWF10x/22mm,Diopter Adjustable			
Zoom Body	0.8~5.0x, Zoom Ratio 1:6	o		
	0.8~6.4x, Zoom Ratio 1:8		o	
	0.8~8.0x, Zoom Ratio 1:10			o
Auxiliary Lens	Plan Achromatic Lens 1.0x			
Working Distance	78mm			
Focusing	Coarse Focusing Konb, Range 105mm	o	o	
	Coaxial Coarse & Fine Focusing Knob, Fine Division 0.015mm, Coarse Stroke 38.4mm/Rotation, Fine Stroke 1.5mm/Rotation, Moving Range 105mm			o
Illumination	Transmit & Reflect LED Illumination, Brightness Adjustable, Working Life 6000 Hours			
A23.1001 Zoom Stereo Microscope		Cata. No.		
Head	Tilting Head 0°-35° Viewing Angle Adjustable for User Change to A23.1005	A53.1003-2		
Eyeiece	WF15x/16mm	A51.1027-1516		
	WF20x/12mm	A51.1027-2012		
	WF30x	A51.1027-30		
Auxiliary Lens	Plan Apochromatic Lens 0.5x	A52.1023-05A		
	Plan Achromatic Lens 1.0x	A52.1023-10		
	Plan Achromatic Lens 2.0x	A52.1023-20		
	Achromatic Lens 0.3x	A52.1023-03		
	Achromatic Lens 0.5x	A52.1023-05		
Illumination	Coaxial Illumination Set: --Cold Light Power Supplier, --Dual Light Guide, --Coaxial Episcopic Illuminator, --Iris Diaphragm	A56.1021		
	Iris Diaphragm	A56.1021-2		
Working Stage	Temperature Control Stage	A54.1004		
	Movable Working Stage	A54.1007		
Stand Click to See	Universal Table Stand	A54.1011		
		A54.1012		
		A54.1013		
		A54.1014		
	Track Stand, Transmit & Reflect Light Source With Coarse & Fine Focusing Optional, Model - F, Like A54.1019-F	A54.1019		
		A54.1020		
		A54.1021		
Fluorescent	EPI-Fluorescent Attachment For Stereo Microscope Change to A23.1004	A5F.1020		
Dark Field	Dark Field Attachment For Stereo Microscope	A5D.1001		
Polarizing	Simple Polarizing Attachment For Stereo Microscope	A5P.1003		
LCD Monitor	8.4" TFT Touch Screen LCD Monitor & Connector Change to A36.1001	A59.1006		
Adapter	Photo Tube, Beamsplitter 1 Port	A55.1003-S1		
	Photo Tube, Beamsplitter 2 Port	A55.1003-S2		
	Photo Adapter, 2.0x, For Nikon or Canon Dlgital SLR Camera	A55.1020-2X		
	C-mount 0.5x	A55.1021-05		
	C-mount 1.0x	A55.1021-10		



A23.1001 Parallel Zoom Stereo Microscope




Parallel Optical System, Module Design to Add Function Easily



Full Range Accessories Optional For A23.1001






A23.1001 Parallel Zoom Stereo Microscope System Size



A5F.1020 EPI-Fluorescent Attachment
For Stereo Microscope



A23.1001
Parallel Zoom Stereo
0.8-8x



A23.1004
Fluorescent Zoom
Stereo
Microscope 0.8x-8x



A54.1004 Temperature Control Stage



A54.1007 Movable Working Stage



Black/White Plate Diameter 10cm



Glass Plate Diameter 10cm Optional



Opto-Edu (Beijing) Co., Ltd.

0086 13911110627 sale@optoedu.com cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China