



Digital Microscope Cameras Digital Camera , HDMI , 1080p Microscope Accessory A59.3511

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A59.3511
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Microscope
- Certification: CE|Rohs
- Product Name: Microscope Accessories
- Resolution: 2.0M
- Sensor: CMOS
- Pixel Size: 2.75*2.75um
- Highlight: compound microscope parts,
parts of the microscope



More Images



Product Description

Digital Microscope Cameras Digital Camera, HDMI, 1080p Microscope Accessory A59.3511

A59.3511 Digital Camera, HDMI, 1080p	
Sensor	1/2.8", 1920*1080, Panasonic Sensor
Pixel Size	2.75*2.75um
White balance	Manual /Auto can save white balance parameters
Digital zoom	10X, and narrow function
Light Control	Manual /Auto
Colour	R/G/B Adjustable
Frozen	Support
OSD	UI(User Interface), mouse
Edge Enhancement	Support
Photo and Video	Support
Crosshairs	Support, many colours available and Thickness Adjustable
Image Contrast	Support
Wide Dynamic	Support
Digital Noise Reduction	Support
	HDMI
Interface	USB2.0
	SD Card Slot
Mouse	Wireless Mouse
System	Windows XP/Vista/ Win7(32 and 64-bit)



A59.3511 HDMI Digital Camera



A59.3511 HDMI Digital Camera



www.cnoec.com
A59.3511 HDMI Digital Camera



A59.3511 HDMI Digital Camera





Lens is Optional



Packing Box



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China