



9" LCD Digital Binocular Head Microscope Accessories A59.3701

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A59.3701
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Highlight: microscope accessory, compound microscope parts



More Images



Product Description

9" LCD Digital Binocular Head Microscope Accessories A59.3701

- 9" LCD Screen
- HDMI Output
- With 4G SD Card
- Having 9 languages for choice.

A59.3701 9" LCD Digital Binocular Head	
Eyepiece	WF10x/18mm
Pixel	5Mega
Photo Resolution	2560X1920
Video Resolution	1280X720/30FPS
LCD Screen	9 HD Inches LCD Screen, Resolution is 1280X 800
Data Output	HDMI, USB (For Win XP),SD Slot
Storage	4G SD Card
Exposure Mode	Auto Exposure
Menu Language	Nine Languages (English,Germany,French,Spanish,Russia,Janpanese,Korean,La tin,Arabic)



A59.3701






9 HD Inches LCD and WF10x/18mm



Side View of A59.3701



Interfaces





Backside View of A59.3701



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China